



Yes, you can:

Tel: 0800 043 0881
Email: salesUK@dhcare.com
Online Spares available at: www.dhcare.co.uk

Power Retail Prescription Form

September 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX20 Modulite

Standard Features

- Max User Weight 136kg (21 Stone)
- Crash tested ISO7176 - 19
- Highly adjustable Modulite Seating System
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths from 410mm - 510mm (16" - 20"), stepless adjustment
- Tension adjustable backrest
- Standard T-Armrests, height and depth adjustable, detachable, with short armpads
- Manual settable Seat Tilt (0° - 15°)
- Manual settable Backrest Recline (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplates
- LiNX DLX-REM110 with standard bracket
- Choice of 5 shroud colours and 5 rim insert colours
- 8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension (Control - Traction - Comfort)
- 6 km/h (4mph) 4 pole motors and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £6,500



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code Description

BASE MODEL*

JBC0000		AVIVA RX20 Modulite	£6,500
JBZ1300		Special AVIVA RX20 Modulite (Extended Lead Time). Please provide details in comments section.	POA

CHASSIS WIDTH*

JBC2810		Narrow (560 mm/22")	◆
UNZ2830		Wide (635 mm/25") (NC 0000)	POA

SEATING TYPE*

JBC2100		Modulite Seat Unit (NC 2400, 2401, 2402, 2403, 4520, 4521, 2410, 2411, 2412, 2328, 2329, 2330, 2331, 4510, 4530, 4540, 4550)	◆
JBC2103		Modulite Support+ (NC 2112, 2104, 2132, 2140, 2142, 2143, 2144, 2147, 2148, 1980, 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985, 2163, 2170, 2171, 2172, 2173, 2174, 2175, 2176, 2177, 2178, 2168, 2169, 2179, 2380, 2381, 2382, 2383, 2385, 2386, 2387, 2389, 2390, 2391, 2393, 2394, 2395, 2151, 2152, 2155, 2156, 2195, 2196, 2197, 2199, 2184, 2186, 2188, 2185, 2187, 2189, 1350, 1427, 1432, 1495, 1496, 2324, 2325, 2327)	○

SEAT FRAME*

JBC2111		Telescopic Seat Frame + Plate with Side Rails. Width 380 mm to 530 mm. Additional 50 mm width by arm rest adjustment. (NC 2156)	◆
JBC2112		Telescopic Seat Frame + Sling Seat with Side Rails. Width 380 mm to 530 mm. +30 mm seat height. Additional 50 mm width by arm rest adjustment. (NC 2156)	○
JBX2102		Telescopic Seat Frame and Plate with Side Rails width 530 to 550 mm (NC 2122, 2123)	£110

SEAT WIDTH*

JBC2121		380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests. (NC 2174, 2175, 2176, 2148)	○
JBC2122		430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.	○
JBC2123		480 mm Seat Plate (NC 2704, 2706, 2148)	○
JBC2124		530 mm Seat Plate (NC 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706, 2106, 2148)	○
JBC2104		550 mm Seat Plate (NC 2132, 2134, 2142, 2143, 2144, 2147, 1980, 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985, 2172, 2173, 2174, 2175, 2176, 2177, 2178, 2168, 2169, 2179, 2380, 2381, 2382, 2383, 2385, 2386, 2387, 2389, 2390, 2391, 2393, 2394, 2395, 2151, 2152, 2155, 2163, 2192, 2193, 2194, 2198, 2199, 2184, 2186, 2188, 2185, 2187, 2189, 1496, 1445, 0556, 0589, 0526, 2701, 2702, 2703, 2704, 2705, 2706, 1404, 1405)	○

BODY FIT KIT

JBC2106		Body Fit Kit - Outrigging Backrest Brackets and Pivot Point (outrig: 1" /25mm per side) (NC 2192, 2194)	○
---------	--	---	---

SEAT DEPTH*

JBC2132		410 mm (NC 2197, 2198, 2701, 2702, 2703, 2704, 2705, 2706)	○
JBC2134		460 mm	○
JBC2136		510 mm	○

SEAT CUSHION*

JBC2142		Flo-tech Lite Cushion (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	£120
JBC2143		Flo-tech Contour Visco; (2-way stretch material) (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	£169
JBC2147		Matrx Libra Cushion, 90 mm thick (Extended lead time)	£410
JBC2144		Matrx Vi Cushion (Mesh fabric) (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983)	£400
JBC2148		Standard Cushion (90mm high) - For 550mm Seat Plate Only (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	○
JBC2149		Support+ Cushion	£270
JBC2140		Without seat cushion (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	○

1st POSITIONING ACCESORIES FOR MATRX LIBRA			
JBC1980		Tapered Wedges (Pair)	£30
JBC1981		Abductor	£23
JBC1982		Lateral Wedge	£20
JBC1983		Anterior Seat Wedges	£27
2nd POSITIONING ACCESORIES FOR MATRX LIBRA			
JBC1980		Tapered Wedges (Pair)	£30
JBC1981		Abductor	£23
JBC1982		Lateral Wedge	£20
JBC1983		Anterior Seat Wedges	£27
POSITIONING ACCESSORIES FOR MATRX VI			
JBC1985		Pelvic Obliquity Kit (Contains 2x positioning accessories)	£54
BACKREST SUPPORT*			
JBC2163		Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. Push Handles included.	(NC 2173, 2052, 2053) ○
JBC2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. Push Handles included.	○
PUSH HANDLES FOR BACKPOSTS			
JBC2108		Backposts without Push Handles	(NC 2219, 2220, 2227, 2228, 2245, 2291, 2292) ○
BACKREST UPHOLSTERY / CUSHION*			
JBC2171		Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.	(NC 2052, 2053) ◆
JBC2179		Flex3 Tension Adjustable Upholstery. Height 540 mm (seat plate to top of upholstery).	(NC 2192, 2193, 2194, 2199, 2155) £155
JBC2172		Flex3 adjustable Backrest Plates (without Backrest Cushion). Flex3 Backrest plate width and height adjustable 560 mm - 790 mm	(NC 2128) £250
JBC2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion. Flex3 Backrest plate width and height adjustable 560 mm - 790 mm	(NC 2128) £400
JBC2400		Support+ Standard Backrest (Spacer fabric)	(NC 2152, 2155, 2156) £330
JBC2401		Support+ Active Strapped Backrest (Spacer fabric)	(NC 2152, 2155, 2156) £400
JBC2402		Support+ Passive Strapped Backrest (Spacer fabric)	(NC 2152, 2155, 2156) £400
JBC2403		Support+ Extra Deep Backrest (Spacer fabric)	(NC 2152, 2155, 2156) £450
JBC2174		Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen	£353
JBC2175		Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	£353
JBC2176		Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen	£494
JBC2177		Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	£558
JBC2178		Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	£546
JBC2168		Matrx MAC Back - Foam Cushion [Mesh fabric]	(NC 1445) £475
JBC2169		Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]	(NC 1445) £575
JBC2170		Without Backrest Upholstery (Back Posts/Back Frame only)	(NC 2324, 2325, 2327, 2151, 2152, 2155) ○

SUPPORT+ BACKREST HEIGHT (Only applicable for Support+ Back options)				
JBC4520		Tall - 570 mm (Measured from Seat Plate)		○
JBC4521		Short - 500 mm (Measured from Seat Plate)		○
MAC BACK BACKREST HEIGHT (Only applicable for MAC Back options)				
JBC2380		405 mm (16 inch) MAC BH	(NC 1445)	○
JBC2381		455 mm (18 inch) MAC BH	(NC 1445)	○
JBC2382		505 mm (20 inch) MAC BH	(NC 1445)	○
JBC2383		560 mm (22 inch) MAC BH	(NC 1445)	○
MAC BACK LATERAL CONTOUR DEPTH LH (Only applicable for MAC Back options)				
JBC2385		90 mm (3.5 inch) MAC Contour Depth LH	(NC 1445)	£100
JBC2386		150 mm (6 inch) MAC Contour Depth LH	(NC 1445)	£100
JBC2387		230 mm (9 inch) MAC Contour Depth LH	(NC 1445)	£100
MAC BACK LATERAL CONTOUR DEPTH RH (Only applicable for MAC Back options)				
JBC2389		90 mm (3.5 inch) MAC Contour Depth RH	(NC 1445)	£100
JBC2390		150 mm (6 inch) MAC Contour Depth RH	(NC 1445)	£100
JBC2391		230 mm (9 inch) MAC Contour Depth RH	(NC 1445)	£100
MAC BACK BACKREST WIDTH (Only applicable for MAC Back options)				
JBC2393		Backrest -25 mm (-1 inch) to seat width	(NC 1445)	○
JBC2394		Backrest matches seat width	(NC 1445)	○
JBC2395		Backrest +25 mm (+1 inch) to seat width	(NC 1445)	○
UPHOLSTERY COLOUR*				
JBC2151		Black Nylon Fabric TR18		◆
JBC2152		Black Multi Stretch PU Coated Wipe Down Fabric TR26		○
JBC2155		Black Spacer Fabric TR36	(NC 2171)	○
JBC2156		Black Nylon (For 550mm SW Only)		○
JBX1002		No Upholstery Colour required - No backrest upholstery required	(NC 2179, 2173, 2174, 2175, 2176, 2177, 2178)	○
ARMRESTS*				
JBC2191		T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm		◆
JBC2181		T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm		£40
JBC2198		Padded Armrest with padded side guard (Water Fall Type)	(NC 2195, 2196, 2197, 2199)	£120
JBC2183		Inrigged Armrests. Reduces SW by 25 mm each side		£30
JBC2193		Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm).	(NC 1445, 2179)	£65
JBC2192		Following Armrest. Height adjustable from 230 - 300 mm	(NC 2179, 2106)	£185
JBC2194		Following Armrest. Height adjustable from 300 - 360 mm	(NC 2179, 2106)	£185
ARMPADS*				
JBC2195		Short Armpad (325 mm)		◆
JBC2196		Long Armpad (375 mm)	(NC 2198)	○
JBC2197		Long Ergonomic Armpad (397 mm)	(NC 2198)	○
JBC2410		Support+ Standard Armpad Short (80 mm x 300 mm)		£100
JBC2411		Support+ Standard Armpad Long (80 mm x 400 mm)		£110
JBC2412		Support+ P-Shaped Armpad (120 mm x 400 mm)		£130
JBC2199		Soft Armpad (REA Azalea Wide Pad)		£95
JBX3000		Please select with Option 2198 - Padded Armrest with padded side guard (Water Fall Type)	(NC 2181, 2183, 2193, 2192)	○

HIP SUPPORTS LH			
JBC2184		Hip Support LH, pad size 140 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
JBC2186		Hip Support LH, pad size 230 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
JBC2188		Hip Support LH, pad size 330 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
HIP SUPPORTS RH			
JBC2185		Hip Support RH, pad size 140 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
JBC2187		Hip Support RH, pad size 230 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
JBC2189		Hip Support RH, pad size 330 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
MOUNTING HARDWARE FOR HIP SUPPORTS			
JBC2164		Mounting hardware for hip supports (6,1" / 155 mm high)	○
JBC2165		Mounting hardware for hip supports (8,1" / 206 mm high)	○
HEADREST / NECKREST (Adequate headrest option must be selected)*			
JBC1350		Neckrest depth adjustable	£80
JBC1432		Headrest sideways adjustable	£225
JBC1427		Headrest with cheek support	£225
JBC2328		Support+ Flex Headrest (Medium). Complete with Mounting Hardware & Stem (_2332).	£250
JBC2329		Support+ Flex Headrest (Large). Complete with Mounting Hardware & Stem (_2332).	£270
JBC2330		Support+ Standard Headrest (Medium). Complete with Mounting Hardware & Stem (_2332).	£230
JBC2331		Support+ Standard Headrest (Large). Complete with Mounting Hardware & Stem (_2332).	£250
JBC2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware (NC 2170)	£325
JBC2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware (NC 2170)	£325
JBC2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware (NC 2170)	£440
JBC1496		Non Invacare Headrest but profile mount fitting Flex3 hard back (NC 2171, 2174, 2175, 2176, 2177, 2178)	○
JBC1495		Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	○
SUPPORT+ ACCESSORIES			
JBC4510		Support+ Wipedown Cushion Cover (PU Coated Type)	£30
JBC4530		Support+ Wipedown Backrest Cover (PU Coated Type)	£30
JBC4540		Support+ Wipedown Armpad Covers (PU Coated Type)	£30
JBC4550		Support+ Wipedown Headrest Cover (PU Coated Type)	£30
LATERAL SUPPORTS			
JBC1445		Lateral Support with extended arm, multi-adjustable pad (NC 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194)	£165

MANUAL / POWERED SEATING OPTIONS*				
JBC1091		Manually settable Seat Tilt 0 - 15°		○
JBC1092		Powered "fixed pivot" Tilt 20°	(NC 2201)	£400
JBC2055		Powered "fixed pivot" Tilt 20° & Sit Smarter	(NC 2201)	£540
JBC1093		Powered Tilt (0° - 30°) with CoG shift	(NC 2201)	£1,100
JBC2057		Powered Tilt (0° - 30°) with CoG shift & Sit Smarter	(NC 2201)	£1,490
JBC2052		Powered Tilt (0° - 50°) with CoG shift	(NC 2201, 2171)	£1,450
JBC2059		Powered Tilt (0° - 50°) with CoG shift & Sit Smarter	(NC 2201, 2171)	£2,000
JBC1094		Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift. 6 point Tie downs must be selected	(NC 2201, 2220, 1910, 1912)	£1,900
JBC2058		Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift & Sit Smarter. 6 point Tie downs must be selected	(NC 2201, 2220, 1910, 1912)	£2,570
JBC2053		Seat Lifter with Powered Tilt (300mm, 45°) with CoG. 6 point Tie downs must be selected	(NC 2201, 2220, 1910, 1912, 2171)	£2,500
JBC2060		Seat Lifter with Powered Tilt (300mm, 45°) with CoG & Sit Smarter. 6 point Tie downs must be selected	(NC 2201, 2220, 1910, 1912, 2171)	£3,370
SEAT TO GROUND HEIGHT*				
JBC2857		435 mm		○
JBC2858		460 mm		○
JBC2859		485 mm		○
BACKREST RECLINER*				
JBC0975		Manual Settable Recliner (in 3.75° steps) 90° - 120°		◆
JBC0977		Powered Recliner 30°		£500
JBC0979		Powered Recliner 30° & Sit Smarter		£670
LEG REST TYPE*				
JBC0556		Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm	(NC 0517, 0512, 0504, 0513, 1244, 1416, 0501)	£370
JBC0589		Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm	(NC 2201, 0512, 0504, 0513, 1244, 1416, 0501)	£1,620
JBC0600		Centre Mounted: LNX Powered Centre Mounted Legrest & Sit Smarter, length compensating. Length 340 - 410 mm	(NC 2201, 0512, 0504, 0513, 1244, 1416, 0501)	£2,190
JBC0526		Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm	(NC 2365, 0559, 0516, 0517, 0501)	◆
JBC0503		Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools	(NC 0559, 0516, 0517, 0554, 0555, 0501)	£120
JBC0533		Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad	(NC 2365, 0559, 0516, 0517, 0554, 0555, 0501)	£235
JBC0531		Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad	(NC 2201, 2365, 0559, 0516, 0517, 0554, 0555, 0501)	£660
JBC0530		Side Mounted: ADE pair & Sit Smarter. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad	(NC 2201, 2365, 0559, 0516, 0517, 0554, 0555, 0501)	£900
JBC0501		Without Leg Rests	(NC 0504, 0513, 2365, 0559, 0516, 0517)	○

FOOTPLATE TYPE*			
JBC2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")	£105
JBC0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	○
JBC0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5)	£135
JBC0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")	£270
JBC0512		Fixed Footplate pair (NC 0533, 0531)	◆
JBC0504		Angle Adjustable Footplate pair (NC 0526)	£75
JBC0513		Angle and depth adjustable Footplate pair	£120
JBC0501		Without footplates for option 0501 without Leg Rests (NC 0503, 0533, 0531)	○
AMPUTEE SUPPORT			
JBC0554		Amputee Support RH	£195
JBC0555		Amputee Support LH	£195
LINX CONTROLS*			
JBC2201		DLX-REM110 (Drive Only) (NC 1092, 1093, 1094, 0977, 0589, 2231, 2232, 2219, 2245, 2270, 1908)	◆
JBC2203		DLX-REM211 (Drive & Actuators) (NC 1437, 2231, 2232, 2270, 1908)	○
JBC2204		DLX-REM216 (Drive, Actuators & Lights) (NC 2232, 2270)	£200
JBC2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets) (NC 2220)	£425
JBC2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick). (NC 2220, 2270, 1709, 1710, 1207, 1208)	£765
REM400 TOGGLE OPTION			
JBC2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL) (NC 2203, 2204, 2206)	£315
CONTROL BRACKETRY*			
JBC1207		Standard Bracket RH (NC 2207, 2265, 1404)	○
JBC1208		Standard Bracket LH (NC 2208, 2266, 1405)	○
JBC1707		Swing Away RH Bracket (NC 2207, 2265, 1404)	£135
JBC1708		Swing Away LH Bracket (NC 2208, 2266, 1405)	£135
JBC1709		Matrx Resolve Swing in/out Bracket RH (NC 2207, 2265, 1404)	£245
JBC1710		Matrx Resolve Swing in/out Bracket LH (NC 2208, 2266, 1405)	£245
JBC2213		Centre Mount Swing Away to RH Nucleus Bracket (NC 2207, 2265, 1415, 2217, 1404, 1405)	£495
JBC2214		Centre Mount Swing Away to LH Nucleus Bracket (NC 2208, 2266, 1339, 2218, 1404, 1405)	£495
JBC2209		REM500 Swing Away RH Multi-Adjustable Bracket (NC 2203, 2204, 2205, 2207, 2265, 1415, 2217, 1404)	○
JBC2210		REM500 Swing Away LH Multi-Adjustable Bracket (NC 2203, 2204, 2205, 2208, 2266, 1339, 2218, 1405)	○
JBC2227		Bracket RH mounted for Attendant (NC 2220, 2128, 2108, 2219, 1221, 1312)	£65
JBC2228		Bracket LH mounted for Attendant (NC 2220, 2128, 2108, 2219, 1222, 1363)	£65
JBZ2239		Individual mounting option per request (Special). Requires additional information/specification.	POA
DUAL CONTROLS			
JBC2220		LiNX DLX-REM050 (Dual Control used as secondary remote) (NC 2273, 2108, 2205, 2206, 1221, 1222, 2227, 2228, 9260)	£395
JBC2219		LiNX DLX-ACU200 (Compact Dual Control) (NC 2108, 1221, 1222, 2227, 2228, 9260)	£395
JBC2245		LiNX DLX-IDC (Intuitive Dual Control) (NC 2108, 2259, 2260, 9260)	£745
JBC2420		mo-vis Scoot Control w/ dedicated bracket included. (NC 2108, 2259, 2260, 2259, 2260, 1221, 1222)	£3,420
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)			
JBC2259		Dual Control Standard Bracket RH (NC 2291, 2245, 1312)	○
JBC2260		Dual Control Standard Bracket LH (NC 2292, 2245, 1363)	○
JBC1221		IDC Dual Control Detachable Bracket RH (NC 2291, 2220, 2219, 2227, 1312)	○
JBC1222		IDC Dual Control Detachable Bracket LH (NC 2292, 2220, 2219, 2228, 1363)	○
JBC9260		Bracket for mo-vis Scoot (NC 2219, 2245)	○

COMPACT REMOTES				
JBC2231		LiNX DLX-CR400 (Compact Remote)	(NC 2203, 2204, 2220, 3300, 3301, 3310, 3311, 3320, 3330)	£825
JBC2232		LiNX DLX-CR400LF (Compact Remote Low Force)	(NC 2203, 2204, 2220, 3300, 3301, 3310, 3311, 3320, 3330)	£1,075
JBC3000		mo-vis Micro LiNX	(NC 2201, 2203, 2204)	£2,960
JBC3001		mo-vis Multi LiNX	(NC 4401, 4402, 2201, 2203, 2204)	£2,640
JBC3110		mo-vis All Round Low Force Omni	(NC 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£2,200
JBC3111		mo-vis All Round Omni	(NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£1,600
JBC3112		mo-vis Heavy Duty Omni	(NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4300, 4301, 4302, 4400, 4401, 4402, 2201, 2203, 2204)	£3,000
COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)				
JBC2207		RH Mount	(NC 1707, 2209, 2213, 2217, 1404)	○
JBC2208		LH Mount	(NC 1708, 2210, 2214, 2218, 1405)	○
JBC2265		Centre Mount Swing Away to RH (Nucleus)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£495
JBC2266		Centre Mount Swing Away to LH (Nucleus)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£495
JBZ2215		Individual mounting option per request (Special). Requires additional information/specification.		POA
JBC3300		Maxx Resolve RH	(NC 4301, 4302)	£270
JBC3301		Maxx Resolve LH	(NC 4301, 4302)	£270
JBZ3310		BodyPoint Trilock Midline RH	(NC 4301, 4302, 0000)	£480
JBZ3311		BodyPoint Trilock Midline LH	(NC 4301, 4302, 0000)	£480
JBZ3320		Tray Mounted	(NC 4301, 4302, 0000)	○
JBZ3330		Individual mounting option per request (Special). Requires additional information/specification.		POA
CHIN CONTROLS				
JBC3000		mo-vis Micro LiNX	(NC 2201, 2203, 2204)	£2,960
JBC3001		mo-vis Multi LiNX	(NC 4401, 4402, 2201, 2203, 2204)	£2,640
JBC3110		mo-vis All Round Low Force Omni (consider option 3810 for switch access)	(NC 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£2,200
JBC3111		mo-vis All Round Omni (consider option 3810 for switch access)	(NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£1,600
CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)				
JBC3400		mo-vis Chin Control Harness (One Size)	(NC 4301, 4302)	£1,000
JBC3410		mo-vis Manual Lateral Swing Away RH Mounted	(NC 4301, 4302)	£1,520
JBC3411		mo-vis Manual Lateral Swing Away LH Mounted	(NC 4301, 4302)	£1,520
JBC3412		mo-vis Manual Vertical Swing Away RH Mounted	(NC 4301, 4302)	£1,520
JBC3413		mo-vis Manual Vertical Swing Away LH Mounted	(NC 4301, 4302)	£1,520
JBC3420		mo-vis Powered Lateral Swing Away RH Mounted	(NC 4301, 4302)	£2,400
JBC3421		mo-vis Powered Lateral Swing Away LH Mounted	(NC 4301, 4302)	£2,400
JBC3422		mo-vis Powered Vertical Swing Away RH Mounted	(NC 4301, 4302)	£2,400
JBC3423		mo-vis Powered Vertical Swing Away LH Mounted	(NC 4301, 4302)	£2,400
JBZ3330		Individual mounting option per request (Special). Requires additional information/specification.		POA
HEAD CONTROLS				
JBC3200		mo-vis Proportional Head Control Omni	(NC 1401, 2894, 2201, 2203, 2204)	£4,500
JBC3221		mo-vis Switched Head Control Omni	(NC 1401, 2894, 2201, 2203, 2204)	£4,500
JBC3222		mo-vis Switched Head Control w/ Sip & Puff Omni	(NC 1401, 2894, 2201, 2203, 2204)	£4,750
FOOT CONTROL				
JBC3010		mo-vis Foot Control LiNX RH Side	(NC 2203, 2204, 2220)	£3,300
JBC3011		mo-vis Foot Control LiNX LH Side	(NC 2203, 2204, 2220)	£3,300
LiNX 10-WAY SWITCH				
JBC2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access		£500

LiNX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)				
JBC2217		Bracket RH, mounted for user	(NC 2207, 2213, 2209, 2265, 1404)	○
JBC2218		Bracket LH, mounted for user	(NC 2208, 2214, 2210, 2266, 1405)	○
JBC2291		Bracket RH, mounted for attendant	(NC 2259, 1221, 2108)	£65
JBC2292		Bracket LH, mounted for attendant	(NC 2260, 1222, 2108)	£65
LiNX USB CHARGER PORT				
JBC2253		USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706)	£160
LiNX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)				
JBC2281		RH Mount		○
JBC2282		LH Mount		○
LiNX INPUT MODULE				
JBC2273		DLX-IN200 Module (Jack port for external power button)		£270
UNZ9249		DLX-IN500 Module (for DB9 Drive input. Only compatible with DLX-REM400 / DLX-REM500)	(NC 2201, 2203, 2204, 2220)	£785
FURTHER LiNX / MO-VIS COMPONENTS / ACCESSORIES				
JBC2302		Not Necessary	(NC 2254, 1711, 1712, 1713, 1714, 1715, 1716)	◆
JBC2300		To be enclosed only. Standard 300 mm Cable Length for mo-vis Switches		○
JBC2305		To be enclosed only. Extended 1500 mm Cable Length for mo-vis Switches		○
JBZ2301		To be connected and programmed (Special)		POA
BUDDY BUTTONS				
JBC1685		Buddy Button. Red / Yellow / Green / Blue / Black provided. Interchangeable colours.		£90
EGG SWITCH				
JBC2254		Egg Switch, black		£90
PIKO BUTTONS				
JBC1711		Red Piko Button		£160
JBC1712		Yellow Piko Button		£160
JBC1713		Green Piko Button		£160
JBC1714		Orange Piko Button		£160
JBC1715		Blue Piko Button		£160
JBC1716		Black Piko Button		£160
MO-VIS TWISTER BASIC				
JBC3500		Black		£140
JBC3501		Red		£140
JBC3502		Yellow		£140
JBC3503		Blue		£140
JBC3504		Green		£140
MO-VIS TWISTER BASIC MOUNTING (Only applicable for mo-vis Twister Basic)				
JBC3520		Mounting Set		£50
JBC3525		Table Mounting Set		£50
MO-VIS TWISTER				
JBC3600		Black		£160
JBC3601		Red		£160
JBC3602		Yellow		£160
JBC3603		Blue		£160
JBC3604		Green		£160

MO-VIS TWISTER ON BENDED TUBE INC. D6 CLAMP			
JBC3700		Black	£220
JBC3701		Red	£220
JBC3702		Yellow	£220
JBC3703		Blue	£220
JBC3704		Green	£220
MO-VIS SATELLITE TWISTER (Only applicable for mo-vis All Round Remotes)			
JBC3800		Black	£200
JBC3801		Red	£200
JBC3802		Yellow	£200
JBC3803		Blue	£200
JBC3804		Green	£200
JBC3810		All Round Chin Control Set (includes All Round Joystick Ball and 2x Satellite Switches)	£420
MO-VIS TWISTER GOOSENECK			
JBC3900		Black	£260
JBC3901		Red	£260
JBC3902		Green	£260
MO-VIS Q2M TWIST & SWING			
JBC4000		Black	£500
JBC4001		Red	£500
JBC4002		Yellow	£500
JBC4003		Blue	£500
JBC4004		Green	£500
MO-VIS TWISTER PRO			
JBC4100		Black	£160
JBC4101		Red	£160
JBC4102		Green	£160
MO-VIS MULTI SWITCH			
JBC4200		Multi Switch	£620
JBC4201		Multi Switch Proximity Sensors	£860
JBC4202		Multi Switch Proximity Set	£840
MO-VIS HANDWARMER			
JBC4300		Hand Warmer	£800
JBZ4301		Hand Warmer 2 Dome Hand Hood	£400
JBZ4302		Hand Warmer 1 Dome Hand Hood	£300
MO-VIS JOYSTICK KNOBS			
JBC4400		Micro Joystick Ball	£60
JBC4401		Micro Joystick Cup 10.1mm	£60
JBC4402		Micro Joystick Cup 15.5mm	£60
JBC4403		Multi Joystick Ball 4mm	£60
JBC4404		Multi Joystick Cup 24mm	£60
JBC4405		All Round Joystick Ball 6.35mm	£60
JBC4406		All Round Joystick Cup 37mm	£60

BODYPOINT JOYSTICK KNOBS			
JBC2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A	£75
JBC2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A	£75
JBC2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£75
JBC2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£75
5:1 SWITCH ADAPTOR			
JBC1681		5:1 Switch Adapter (9-pin sub D to 5 jack ports)	£120
SHROUD COLOUR*			
JBC2871		Mid range - Black/Anthracite Shroud	○
JBC2872		Mid range - Silver Shroud	○
JBC2873		Mid range - Bengal Red Shroud	○
JBC2874		Mid range - Mercury Blue Shroud	○
JBC2875		Mid range - Midnight Purple Shroud	○
JBZ0299		Premium - Individual Shroud Colour. RAL number required Extended lead lime: + 15 days.	£300
RIM COLOUR*			
JBC2882		Mid range - Black Rim Inserts	○
JBC2883		Mid range - Silver Rim Inserts	○
JBC2884		Mid range - Bengal Red Rim Inserts	○
JBC2885		Mid range - Mercury Blue Rim Inserts	○
JBC2886		Mid range - Midnight Purple Rim Inserts	○
JBZ2887		Premium - Individual Rim Insert Colour. RAL number required Extended lead lime: + 15 days. Price Included with Premium Shroud Kit	○
FORKS*			
JBC2631		Single sided fork	£250
JBC2632		Double sided fork Without castor rim insert	◆
CASTORS*			
JBC0835		9" (225 mm) x 2.80 (70 mm) Castor, pneumatic	○
JBC0836		9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	○
JBC0838		9" (225 mm) x 2.80 (70 mm) Castor, puncture protected Gel-sealed	○
JBC0828		8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre	◆
CASTOR LOCKS			
JBC1490		Castor Lock	£195
DRIVE WHEELS*			
JBC0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	◆
JBC0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	○
C.T.C SUSPENSION*			
JBC2870		CTC Suspension Control - Traction - Comfort	◆
SUSPENSION SETTING*			
JBC1910		Soft (NC 1094)	◆
JBC1912		Medium (NC 1094)	○
JBC1911		Firm	○
JBC1915		Very firm	○

MOTOR TYPE*			
JBC1294		Indoor setting, reduced speed approx. 3 km/h	○
JBC1282		6 km/h	◆
JBC1908		10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 2203, 1437, 1254) £515
BATTERY TYPE*			
JBC1259		60 Ah Gel	◆
G TRAC			
JBC1215		G-Trac Module	£400
CHARGER*			
JBC1343		Opti Charge Charger 8A (IP44)	◆
JBC1313		Opti Charge Charger 10A (IP44)	£155
ANTI-TIPPERS*			
JBC1205		Anti-tipper	◆
LIGHTS			
JBC1437		LED Light system with indicators 30 Lux headlights, included with 10km/h motors	(NC 2203, 1908) £450
JBC1254		Without Light system	○
KERB RIDER			
JBC1244		Kerb Rider with park position feature	(NC 0556, 0589, 1294) £215
TRANSPORTATION KIT*			
JBC1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.	£580
JBC1493		4-Point tie-down (For wheelchair transportation without user)	○
JBC1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. Not compatible with Seat lifter.	◆
JBC1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory.	(NC 1091, 1092, 1093) ○
EASY ACCESS FRONT TIE DOWNS			
JBC2894		Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down. For transport worthy versions (1494/1401) only	£55
FIXATION BELT*			
JBC1329		Lap belt with Auto Lock Standard mount on backrest	○
JBC2113		Lap belt with Auto Lock Flexible mount in side rails	◆
JBC2701		BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £95
JBC2702		BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £95
JBC2703		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £90
JBC2704		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £90
JBC2705		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £90
JBC2706		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £90

REAR VIEW MIRROR			
JBC1415		Rear view mirror RH (NC 2217, 2209, 2213, 2265)	£45
JBC1339		Rear view mirror LH (NC 2218, 2210, 2214, 2266)	£45
BAG HOOKS			
JBC2128		Bag Hooks (Max load up to 6kg = impact on max user weight) (NC 2172, 2173, 2227, 2228, 2291, 2292)	£80
CANE HOLDER			
JBC1312		Cane holder RH (NC 2227, 2259, 1221)	£45
JBC1363		Cane holder LH (NC 2228, 2260, 1222)	£45
WHEEL FENDERS			
JBC1416		Fender for Castors	£80
JBC1417		Fender for Drive wheels	£160
TRAY			
JBC1404		Tray, Swing Away RH (NC 1707, 2209, 2265, 2266, 2207, 2213, 2214, 2217)	£200
JBC1405		Tray, Swing Away LH (NC 1708, 2210, 2265, 2266, 2208, 2213, 2214, 2218)	£200

TOTAL PRICE			
-------------	--	--	--

Non Conforming Options:	
JBC0000	AVIVA RX20 Modulite
JBC2112	Telescopic Seat Frame + Sling Seat with Side Rails. Width 380 mm to 530 mm. +30 mm seat height. Additional 50 mm width by arm rest adjustment.
JBC2121	380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.
JBC2122	430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.
JBC2123	480 mm Seat Plate
JBC2124	530 mm Seat Plate
JBC2104	550 mm Seat Plate
JBC2106	Body Fit Kit - Outrigging Backrest Brackets and Pivot Point (outrig: 1" /25mm per side)
JBC2132	410 mm
JBC2134	460 mm
JBC2142	Flo-tech Lite Cushion
JBC2143	Flo-tech Contour Visco; (2-way stretch material)
JBC2147	Matrx Libra Cushion, 90 mm thick (Extended lead time)
JBC2144	Matrx Vi Cushion (Mesh fabric)
JBC2148	Standard Cushion (90mm high) - For 550mm Seat Plate Only
JBC2140	Without seat cushion
JBC1980	Tapered Wedges (Pair)
JBC1981	Abductor
JBC1982	Lateral Wedge
JBC1983	Anterior Seat Wedges
JBC1980	Tapered Wedges (Pair)
JBC1981	Abductor
JBC1982	Lateral Wedge
JBC1983	Anterior Seat Wedges
JBC1985	Pelvic Obliquity Kit (Contains 2x positioning accessories)
JBC2163	Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. Push Handles included.
JBC2108	Backposts without Push Handles
JBC2171	Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.
JBC2179	Flex3 Tension Adjustable Upholstery. Height 540 mm (seat plate to top of upholstery).
JBC2172	Flex3 adjustable Backrest Plates (without Backrest Cushion). Flex3 Backrest plate width and height adjustable 560 mm - 790 mm
JBC2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion. Flex3 Backrest plate width and height adjustable 560 mm - 790 mm
JBC2400	Support+ Standard Backrest (Spacer fabric)
JBC2401	Support+ Active Strapped Backrest (Spacer fabric)
JBC2402	Support+ Passive Strapped Backrest (Spacer fabric)
JBC2403	Support+ Extra Deep Backrest (Spacer fabric)
JBC2174	Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen
JBC2175	Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBC2176	Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen

JBC2177	Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBC2178	Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBC2168	Matrx MAC Back - Foam Cushion [Mesh fabric]
JBC2169	Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]
JBC2170	Without Backrest Upholstery (Back Posts/Back Frame only)
JBC4520	Tall - 570 mm (Measured from Seat Plate)
JBC4521	Short - 500 mm (Measured from Seat Plate)
JBC2380	405 mm (16 inch) MAC BH
JBC2381	455 mm (18 inch) MAC BH
JBC2382	505 mm (20 inch) MAC BH
JBC2383	560 mm (22 inch) MAC BH
JBC2385	90 mm (3.5 inch) MAC Contour Depth LH
JBC2386	150 mm (6 inch) MAC Contour Depth LH
JBC2387	230 mm (9 inch) MAC Contour Depth LH
JBC2389	90 mm (3.5 inch) MAC Contour Depth RH
JBC2390	150 mm (6 inch) MAC Contour Depth RH
JBC2391	230 mm (9 inch) MAC Contour Depth RH
JBC2393	Backrest -25 mm (-1 inch) to seat width
JBC2394	Backrest matches seat width
JBC2395	Backrest +25 mm (+1 inch) to seat width
JBC2151	Black Nylon Fabric TR18
JBC2152	Black Multi Stretch PU Coated Wipe Down Fabric TR26
JBC2155	Black Spacer Fabric TR36
JBC2156	Black Nylon (For 550mm SW Only)
JBC2181	T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm
JBC2198	Padded Armrest with padded side guard (Water Fall Type)
JBC2183	Inrigged Armrests. Reduces SW by 25 mm each side
JBC2193	Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm).
JBC2192	Following Armrest. Height adjustable from 230 - 300 mm
JBC2194	Following Armrest. Height adjustable from 300 - 360 mm
JBC2195	Short Armpad (325 mm)
JBC2196	Long Armpad (375 mm)
JBC2197	Long Ergonomic Armpad (397 mm)
JBC2410	Support+ Standard Armpad Short (80 mm x 300 mm)
JBC2411	Support+ Standard Armpad Long (80 mm x 400 mm)
JBC2412	Support+ P-Shaped Armpad (120 mm x 400 mm)
JBC2199	Soft Armpad (REA Azalea Wide Pad)
JBC2184	Hip Support LH, pad size 140 x 100 mm
JBC2186	Hip Support LH, pad size 230 x 100 mm
JBC2188	Hip Support LH, pad size 330 x 100 mm
JBC2185	Hip Support RH, pad size 140 x 100 mm
JBC2187	Hip Support RH, pad size 230 x 100 mm

JBC2189	Hip Support RH, pad size 330 x 100 mm
JBC1350	Neckrest depth adjustable
JBC1432	Headrest sideways adjustable
JBC1427	Headrest with cheek support
JBC2328	Support+ Flex Headrest (Medium). Complete with Mounting Hardware & Stem (_2332).
JBC2329	Support+ Flex Headrest (Large). Complete with Mounting Hardware & Stem (_2332).
JBC2330	Support+ Standard Headrest (Medium). Complete with Mounting Hardware & Stem (_2332).
JBC2331	Support+ Standard Headrest (Large). Complete with Mounting Hardware & Stem (_2332).
JBC2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware
JBC2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware
JBC2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware
JBC1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
JBC1495	Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts
JBC4510	Support+ Wipedown Cushion Cover (PU Coated Type)
JBC4530	Support+ Wipedown Backrest Cover (PU Coated Type)
JBC4540	Support+ Wipedown Armpad Covers (PU Coated Type)
JBC4550	Support+ Wipedown Headrest Cover (PU Coated Type)
JBC1445	Lateral Support with extended arm, multi-adjustable pad
JBC1091	Manually settable Seat Tilt 0 - 15°
JBC1092	Powered "fixed pivot" Tilt 20°
JBC1093	Powered Tilt (0° - 30°) with CoG shift
JBC2052	Powered Tilt (0° - 50°) with CoG shift
JBC1094	Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift. 6 point Tie downs must be selected
JBC2053	Seat Lifter with Powered Tilt (300mm, 45°) with CoG. 6 point Tie downs must be selected
JBC0977	Powered Recliner 30°
JBC0556	Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm
JBC0589	Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm
JBC0526	Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm
JBC0503	Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools
JBC0533	Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad
JBC0531	Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad
JBC0501	Without Leg Rests
JBC2365	Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")
JBC0559	Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")
JBC0516	Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5")
JBC0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
JBC0512	Fixed Footplate pair
JBC0504	Angle Adjustable Footplate pair
JBC0513	Angle and depth adjustable Footplate pair
JBC0501	Without footplates for option 0501 without Leg Rests
JBC0554	Amputee Support RH

JBC0555	Amputee Support LH
JBC2201	DLX-REM110 (Drive Only)
JBC2203	DLX-REM211 (Drive & Actuators)
JBC2204	DLX-REM216 (Drive, Actuators & Lights)
JBC2205	DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)
JBC2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
JBC2270	Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)
JBC1207	Standard Bracket RH
JBC1208	Standard Bracket LH
JBC1707	Swing Away RH Bracket
JBC1708	Swing Away LH Bracket
JBC1709	Matrx Resolve Swing in/out Bracket RH
JBC1710	Matrx Resolve Swing in/out Bracket LH
JBC2213	Centre Mount Swing Away to RH Nucleus Bracket
JBC2214	Centre Mount Swing Away to LH Nucleus Bracket
JBC2209	REM500 Swing Away RH Multi-Adjustable Bracket
JBC2210	REM500 Swing Away LH Multi-Adjustable Bracket
JBC2227	Bracket RH mounted for Attendant
JBC2228	Bracket LH mounted for Attendant
JBC2220	LiNX DLX-REM050 (Dual Control used as secondary remote)
JBC2219	LiNX DLX-ACU200 (Compact Dual Control)
JBC2245	LiNX DLX-IDC (Intuitive Dual Control)
JBC2259	Dual Control Standard Bracket RH
JBC2260	Dual Control Standard Bracket LH
JBC1221	IDC Dual Control Detachable Bracket RH
JBC1222	IDC Dual Control Detachable Bracket LH
JBC9260	Bracket for mo-vis Scoot
JBC2231	LiNX DLX-CR400 (Compact Remote)
JBC2232	LiNX DLX-CR400LF (Compact Remote Low Force)
JBC2207	RH Mount
JBC2208	LH Mount
JBC2265	Centre Mount Swing Away to RH (Nucleus)
JBC2266	Centre Mount Swing Away to LH (Nucleus)
JBC3300	Maxx Resolve RH
JBC3301	Maxx Resolve LH
JBZ3310	BodyPoint Trilock Midline RH
JBZ3311	BodyPoint Trilock Midline LH
JBZ3320	Tray Mounted
JBZ3330	Individual mounting option per request (Special). Requires additional information/specification.
JBC3410	mo-vis Manual Lateral Swing Away RH Mounted
JBC3411	mo-vis Manual Lateral Swing Away LH Mounted
JBC3412	mo-vis Manual Vertical Swing Away RH Mounted
JBC3413	mo-vis Manual Vertical Swing Away LH Mounted
JBC3420	mo-vis Powered Lateral Swing Away RH Mounted

JBC3421	mo-vis Powered Lateral Swing Away LH Mounted
JBC3422	mo-vis Powered Vertical Swing Away RH Mounted
JBC3423	mo-vis Powered Vertical Swing Away LH Mounted
JBZ3330	Individual mounting option per request (Special). Requires additional information/specification.
JBC3200	mo-vis Proportional Head Control Omni
JBC3221	mo-vis Switched Head Control Omni
JBC3222	mo-vis Switched Head Control w/ Sip & Puff Omni
JBC2217	Bracket RH, mounted for user
JBC2218	Bracket LH, mounted for user
JBC2291	Bracket RH, mounted for attendant
JBC2292	Bracket LH, mounted for attendant
JBC2253	USB Charger Port DLX-USB02-A
JBC2273	DLX-IN200 Module (Jack port for external power button)
JBC2254	Egg Switch, black
JBC1711	Red Piko Button
JBC1712	Yellow Piko Button
JBC1713	Green Piko Button
JBC1714	Orange Piko Button
JBC1715	Blue Piko Button
JBC1716	Black Piko Button
JBC4300	Hand Warmer
JBZ4301	Hand Warmer 2 Dome Hand Hood
JBZ4302	Hand Warmer 1 Dome Hand Hood
JBC4400	Micro Joystick Ball
JBC4401	Micro Joystick Cup 10.1mm
JBC4402	Micro Joystick Cup 15.5mm
JBC1910	Soft
JBC1912	Medium
JBC1294	Indoor setting, reduced speed approx. 3 km/h
JBC1908	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBC1437	LED Light system with indicators 30 Lux headlights, included with 10km/h motors
JBC1254	Without Light system
JBC1244	Kerb Rider with park position feature
JBC1401	6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory.
JBC2894	Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down. For transport worthy versions (1494/1401) only
JBC2701	BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm
JBC2702	BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm
JBC2703	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm
JBC2704	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm
JBC2705	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm
JBC2706	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm

JBC1415	Rear view mirror RH
JBC1339	Rear view mirror LH
JBC2128	Bag Hooks (Max load up to 6kg = impact on max user weight)
JBC1312	Cane holder RH
JBC1363	Cane holder LH
JBC1416	Fender for Castors
JBC1404	Tray, Swing Away RH
JBC1405	Tray, Swing Away LH



For a more complex customer requirement,
please email: salesUK@dhcare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):