



Yes, you can:

Tel: 0800 043 0881
Email: salesUK@dhcare.com
Online Spares available at: www.dhcare.co.uk

Power Retail Prescription Form

September 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX40 HD Modulite

Standard Features

- Max User Weight 160kg (25 Stone)
- Crash tested ISO7176 - 19
- Highly adjustable Modulite Seating System
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths from 410mm - 510mm (16" - 20"), stepless adjustment
- Tension adjustable backrest
- Standard T-Armrests, height and depth adjustable, detachable, with short armpads
- Manual settable Seat Tilt (0° - 15°)
- Manual settable Backrest Recline (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplates
- LiNX DLX-REM110 with standard bracket
- Choice of 5 shroud colours and 5 rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension (Control - Traction - Comfort)
- 6 km/h (4mph) 4 pole motors and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £7,250



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
JBH0002	AVIVA RX40 HD Modulite (NC 2193)	£7,250
JBZ1300	Special AVIVA RX40 HD Modulite (Extended Lead Time). Please provide details in comments section.	POA
CHASSIS WIDTH*		
JBH2830	Wide (635 mm/25")	☒

SEATING TYPE*				
JBH2100		Modulite Seat Unit	(NC 2400, 2401, 2402, 2403, 4520, 4521, 2410, 2411, 2412, 2328, 2329, 2330, 2331, 4510, 4530, 4540, 4550)	○
JBH2103		Modulite Support+	(NC 2112, 2104, 2132, 2140, 2142, 2144, 2147, 2148, 1980, 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985, 2163, 2170, 2171, 2172, 2173, 2174, 2175, 2176, 2177, 2178, 2168, 2169, 2179, 2380, 2381, 2382, 2383, 2385, 2386, 2387, 2389, 2390, 2391, 2393, 2394, 2395, 2151, 2152, 2155, 2156, 2195, 2196, 2197, 2199, 2184, 2186, 2188, 2185, 2187, 2189, 1350, 1427, 1432, 1495, 1496, 2324, 2325, 2327)	○
SEAT FRAME*				
JBH2111		Telescopic Seat Frame + Plate with Side Rails. Width 380 mm to 530 mm. Additional 50 mm width by arm rest adjustment.	(NC 2156)	○
JBH2112		Telescopic Seat Frame + Sling Seat with Side Rails. Width 380 mm to 530 mm. +30 mm seat height. Additional 50 mm width by arm rest adjustment.	(NC 2156)	○
JBX2102		Telescopic Seat Frame and Plate with Side Rails width 530 to 550 mm	(NC 2122, 2123)	£110
SEAT WIDTH*				
JBH2121		380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.	(NC 2174, 2175, 2176, 2148)	○
JBH2122		430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.		○
JBH2123		480 mm Seat Plate	(NC 2702, 2704, 2706, 2148)	○
JBH2124		530 mm Seat Plate	(NC 2174, 2175, 2176, 2177, 2178, 2702, 2704, 2706, 2106, 2148)	○
JBH2104		550 mm Seat Plate	(NC 2132, 2134, 2142, 2144, 2147, 1980, 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985, 2172, 2173, 2174, 2175, 2176, 2177, 2178, 2168, 2169, 2179, 2380, 2381, 2382, 2383, 2385, 2386, 2387, 2389, 2390, 2391, 2393, 2394, 2395, 2151, 2152, 2155, 2192, 2193, 2194, 2198, 2199, 2184, 2186, 2188, 2185, 2187, 2189, 1496, 1445, 0556, 0589, 2701, 2702, 2703, 2704, 2705, 2706, 1404, 1405)	○
BODY FIT KIT				
JBH2106		Body Fit Kit - Outrigging Backrest Brackets and Pivot Point (outrig: 1" /25mm per side)		○
SEAT DEPTH*				
JBH2132		410 mm	(NC 2196, 2197, 2198)	○
JBH2134		460 mm	(NC 2148)	○
JBH2136		510 mm		○
SEAT CUSHION*				
JBH2142		Flo-tech Lite Cushion	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	£120
JBH2147		Matrx Libra Cushion, 90 mm thick (Extended lead time)		£410
JBH2144		Matrx Vi Cushion (Mesh fabric)	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983)	£400
JBH2148		Standard Cushion (90mm high) - For 550mm Seat Plate Only	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	○
JBH2149		Support+ Cushion		£270
JBH2140		Without seat cushion	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	○
1st POSITIONING ACCESORIES FOR MATRX LIBRA				
JBH1980		Tapered Wedges (Pair)		£30
JBH1981		Abductor		£23
JBH1982		Lateral Wedge		£20
JBH1983		Anterior Seat Wedges		£27
2nd POSITIONING ACCESORIES FOR MATRX LIBRA				
JBH1980		Tapered Wedges (Pair)		£30
JBH1981		Abductor		£23
JBH1982		Lateral Wedge		£20
JBH1983		Anterior Seat Wedges		£27
POSITIONING ACCESSORIES FOR MATRX Vi				
JBH1985		Pelvic Obliquity Kit (Contains 2x positioning accessories)		£54

BACKREST SUPPORT*				
JBH2163		Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. Push Handles included.	(NC 2172, 2173)	○
JBH2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. Push Handles included.		◆
PUSH HANDLES FOR BACKPOSTS				
JBH2108		Backposts without Push Handles	(NC 2219, 2227, 2228, 2245, 2291, 2292)	○
BACKREST UPHOLSTERY / CUSHION*				
JBH2171		Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.	(NC 2155)	◆
JBH2179		Flex3 Tension Adjustable Upholstery. Height 540 mm (seat plate to top of upholstery).	(NC 2151, 2192, 2193, 2194, 2199, 2155)	£155
JBH2172		Flex3 adjustable Backrest Plates (without Backrest Cushion). Flex3 Backrest plate width and height adjustable 560 mm - 790 mm	(NC 2163, 2128)	£250
JBH2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion. Flex3 Backrest plate width and height adjustable 560 mm - 790 mm	(NC 2163, 2128)	£400
JBH2400		Support+ Standard Backrest (Spacer fabric)	(NC 2152, 2155, 2156)	£330
JBH2401		Support+ Active Strapped Backrest (Spacer fabric)	(NC 2152, 2155, 2156)	£400
JBH2402		Support+ Passive Strapped Backrest (Spacer fabric)	(NC 2152, 2155, 2156)	£400
JBH2403		Support+ Extra Deep Backrest (Spacer fabric)	(NC 2152, 2155, 2156)	£450
JBH2174		Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 2121)	£353
JBH2175		Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 2121)	£353
JBH2176		Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 2121)	£494
JBH2177		Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen		£558
JBH2178		Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen		£546
JBH2168		Matrx MAC Back - Foam Cushion [Mesh fabric]	(NC 1445)	£475
JBH2169		Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]	(NC 1445)	£575
JBH2170		Without Backrest Upholstery (Back Posts/Back Frame only)	(NC 2324, 2325, 2327, 2151, 2152, 2155)	○
SUPPORT+ BACKREST HEIGHT (Only applicable for Support+ Back options)				
JBH4520		Tall - 570 mm (Measured from Seat Plate)		○
JBH4521		Short - 500 mm (Measured from Seat Plate)		○
MAC BACK BACKREST HEIGHT (Only applicable for MAC Back options)				
JBH2380		405 mm (16 inch) MAC BH	(NC 1445)	○
JBH2381		455 mm (18 inch) MAC BH	(NC 1445)	○
JBH2382		505 mm (20 inch) MAC BH	(NC 1445)	○
JBH2383		560 mm (22 inch) MAC BH	(NC 1445)	○
MAC BACK LATERAL CONTOUR DEPTH LH (Only applicable for MAC Back options)				
JBH2385		90 mm (3.5 inch) MAC Contour Depth LH	(NC 1445)	£100
JBH2386		150 mm (6 inch) MAC Contour Depth LH	(NC 1445)	£100
JBH2387		230 mm (9 inch) MAC Contour Depth LH	(NC 1445)	£100

MAC BACK LATERAL CONTOUR DEPTH RH (Only applicable for MAC Back options)				
JBH2389		90 mm (3.5 inch) MAC Contour Depth RH	(NC 1445)	£100
JBH2390		150 mm (6 inch) MAC Contour Depth RH	(NC 1445)	£100
JBH2391		230 mm (9 inch) MAC Contour Depth RH	(NC 1445)	£100
MAC BACK BACKREST WIDTH (Only applicable for MAC Back options)				
JBH2393		Backrest -25 mm (-1 inch) to seat width	(NC 1445)	○
JBH2394		Backrest matches seat width	(NC 1445)	○
JBH2395		Backrest +25 mm (+1 inch) to seat width	(NC 1445)	○
UPHOLSTERY COLOUR*				
JBH2151		Black Nylon Fabric TR18		◆
JBH2152		Black Multi Stretch PU Coated Wipe Down Fabric TR26	(NC 2171)	○
JBH2155		Black Spacer Fabric TR36	(NC 2171)	○
JBH2156		Black Nylon (For 550mm SW Only)		○
JBX1002		No Upholstery Colour required - No backrest upholstery required	(NC 2179, 2172, 2173, 2174, 2175, 2176, 2177, 2178)	○
ARMRESTS*				
JBH2191		T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm		◆
JBH2181		T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm		£40
JBH2198		Padded Armrest with padded side guard (Water Fall Type)	(NC 2132, 2195, 2196, 2197, 2199)	£120
JBH2193		Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm).	(NC 1445, 2179)	£65
JBH2192		Following Armrest. Height adjustable from 230 - 300 mm	(NC 1445, 2179)	£185
JBH2194		Following Armrest. Height adjustable from 300 - 360 mm	(NC 1445, 2179)	£185
ARMPADS*				
JBH2195		Short Armpad (325 mm)	(NC 2198)	◆
JBH2196		Long Armpad (375 mm)	(NC 2132, 2198)	○
JBH2197		Long Ergonomic Armpad (397 mm)	(NC 2132, 2198)	○
JBH2410		Support+ Standard Armpad Short (80 mm x 300 mm)		£100
JBH2411		Support+ Standard Armpad Long (80 mm x 400 mm)		£110
JBH2412		Support+ P-Shaped Armpad (120 mm x 400 mm)		£130
JBH2199		Soft Armpad (REA Azalea Wide Pad)	(NC 2198)	£95
JBX3000		Please select with Option 2198 - Padded Armrest with padded side guard (Water Fall Type)	(NC 2181, 2193, 2192)	○
HIP SUPPORTS LH				
JBH2184		Hip Support LH, pad size 140 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706)	£250
JBH2186		Hip Support LH, pad size 230 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706)	£250
JBH2188		Hip Support LH, pad size 330 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706)	£250
HIP SUPPORTS RH				
JBH2185		Hip Support RH, pad size 140 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706)	£250
JBH2187		Hip Support RH, pad size 230 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706)	£250
JBH2189		Hip Support RH, pad size 330 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706)	£250
MOUNTING HARDWARE FOR HIP SUPPORTS				
JBH2164		Mounting hardware for hip supports (6,1" / 155 mm high)		○
JBH2165		Mounting hardware for hip supports (8,1" / 206 mm high)		○

HEADREST / NECKREST (Adequate headrest option must be selected)*			
JBH1350		Neckrest depth adjustable	£80
JBH1432		Headrest sideways adjustable	£225
JBH1427		Headrest with cheek support	£225
JBH2328		Support+ Flex Headrest (Medium). Complete with Mounting Hardware & Stem (_2332).	£250
JBH2329		Support+ Flex Headrest (Large). Complete with Mounting Hardware & Stem (_2332).	£270
JBH2330		Support+ Standard Headrest (Medium). Complete with Mounting Hardware & Stem (_2332).	£230
JBH2331		Support+ Standard Headrest (Large). Complete with Mounting Hardware & Stem (_2332).	£250
JBH2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware (NC 2170)	£325
JBH2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware (NC 2170)	£325
JBH2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware (NC 2170)	£440
JBH1496		Non Invacare Headrest but profile mount fitting Flex3 hard back (NC 2170, 2171, 2174, 2175, 2176, 2177, 2178)	○
JBH1495		Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	○
SUPPORT+ ACCESSORIES			
JBH4510		Support+ Wipedown Cushion Cover (PU Coated Type)	£30
JBH4530		Support+ Wipedown Backrest Cover (PU Coated Type)	£30
JBH4540		Support+ Wipedown Armpad Covers (PU Coated Type)	£30
JBH4550		Support+ Wipedown Headrest Cover (PU Coated Type)	£30
LATERAL SUPPORTS			
JBH1445		Lateral Support with extended arm, multi-adjustable pad (NC 2170, 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194)	£165
MANUAL / POWERED SEATING OPTIONS*			
JBH1091		Manually settable Seat Tilt 0 - 15° (NC 1401)	○
JBH1092		Powered "fixed pivot" Tilt 20° (NC 2201)	£400
JBH2055		Powered "fixed pivot" Tilt 20° & Sit Smarter (NC 2201)	£540
JBH1093		Powered Tilt (0° - 30°) with CoG shift (NC 2201)	£1,100
JBH2057		Powered Tilt (0° - 30°) with CoG shift & Sit Smarter (NC 2201)	£1,490
JBH1094		Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift. 6 point Tie downs must be selected (NC 2201)	£1,900
JBH2058		Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift & Sit Smarter. 6 point Tie downs must be selected (NC 2201)	£2,570
SEAT TO GROUND HEIGHT*			
JBH2857		435 mm	○
JBH2858		460 mm	○
JBH2859		485 mm	○
BACKREST RECLINER*			
JBH0975		Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
JBH0977		Powered Recliner 30°	£500
JBH0979		Powered Recliner 30° & Sit Smarter	£670

LEG REST TYPE*			
JBH0556		Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm	(NC 0517, 0512, 0504, 0513, 1244, 1416) £370
JBH0589		Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm	(NC 2201, 0512, 0504, 0513, 1244, 1416) £1,620
JBH0600		Centre Mounted: LNX Powered Centre Mounted Legrest & Sit Smarter, length compensating. Length 340 - 410 mm	(NC 2201, 0512, 0504, 0513, 1244, 1416, 0501) £2,190
JBH0503		Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools	(NC 0559, 0516, 0517, 0501) ◆
JBH0533		Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad	(NC 2365, 0559, 0516, 0517, 0501) £235
JBH0531		Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad	(NC 2201, 2365, 0559, 0516, 0517, 0501) £660
JBH0530		Side Mounted: ADE pair & Sit Smarter. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad	(NC 2201, 2365, 0559, 0516, 0517, 0501) £900
JBH0501		Without Leg Rests	(NC 0504, 0513, 2365, 0559, 0516, 0517) ○
FOOTPLATE TYPE*			
JBH2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")	£105
JBH0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	○
JBH0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5")	£135
JBH0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")	£270
JBH0512		Fixed Footplate pair	(NC 0533, 0531) ◆
JBH0504		Angle Adjustable Footplate pair	£75
JBH0513		Angle and depth adjustable Footplate pair	£120
JBH0501		Without footplates for option 0501 without Leg Rests	(NC 0503, 0533, 0531) ○
LINX CONTROLS*			
JBH2201		DLX-REM110 (Drive Only)	(NC 1092, 1093, 1094, 0977, 0589, 2231, 2232, 2219, 2245, 2270, 1908) ◆
JBH2203		DLX-REM211 (Drive & Actuators)	(NC 1437, 2231, 2232, 2270, 1908, 1909) ○
JBH2204		DLX-REM216 (Drive, Actuators & Lights)	(NC 2232, 2270) £200
JBH2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)	£425
JBH2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).	(NC 2270, 1709, 1710) £765
REM400 TOGGLE OPTION			
JBH2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)	(NC 2203, 2204, 2206) £315

CONTROL BRACKETRY*				
JBH1207		Standard Bracket RH	(NC 2208, 2266, 1405)	○
JBH1208		Standard Bracket LH	(NC 2207, 2265, 1404)	○
JBH1707		Swing Away RH Bracket	(NC 2208, 2266, 1405)	£135
JBH1708		Swing Away LH Bracket	(NC 2207, 2265, 1404)	£135
JBH1709		Matrx Resolve Swing in/out Bracket RH	(NC 2208, 2266, 1405)	£245
JBH1710		Matrx Resolve Swing in/out Bracket LH	(NC 2207, 2265, 1415, 2217, 1404, 1405)	£245
JBH2213		Centre Mount Swing Away to RH Nucleus Bracket	(NC 2208, 2266, 1339, 2218, 1404, 1405)	£495
JBH2214		Centre Mount Swing Away to LH Nucleus Bracket	(NC 2203, 2204, 2205, 2207, 2265, 1415, 2217, 1404)	£495
JBH2209		REM500 Swing Away RH Multi-Adjustable Bracket	(NC 2203, 2204, 2205, 2208, 2266, 1339, 2218, 1405)	○
JBH2210		REM500 Swing Away LH Multi-Adjustable Bracket	(NC 2128, 2108, 2219, 1221, 1312)	○
JBH2227		Bracket RH mounted for Attendant	(NC 2128, 2108, 2219, 1222, 1363)	£65
JBH2228		Bracket LH mounted for Attendant		£65
JBZ2239		Individual mounting option per request (Special). Requires additional information/specification.	(NC 2273, 2108, 2205, 2206, 1221, 1222, 2227, 2228)	POA
DUAL CONTROLS				
JBH2219		LiNX DLX-ACU200 (Compact Dual Control)	(NC 2108, 1221, 1222, 2227, 2228, 9260)	£395
JBH2245		LiNX DLX-IDC (Intuitive Dual Control)	(NC 2108, 2259, 2260, 9260)	£745
JBH2420		mo-vis Scoot Control w/ dedicated bracket included.	(NC 2108, 2259, 2260, 2259, 2260, 1221, 1222)	£3,420
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)				
JBH2259		Bracket RH, Attendant	(NC 2291, 2245, 1312)	○
JBH2260		Bracket LH, Attendant	(NC 2292, 2245, 1363)	○
JBH1221		Bracket RH, detachable	(NC 2291, 2219, 2227, 1312)	○
JBH1222		Bracket LH, detachable	(NC 2292, 2219, 2228, 1363)	○
JBH9260		Bracket for mo-vis Scoot	(NC 2219, 2245)	○
COMPACT REMOTES				
JBH2231		LiNX DLX-CR400 (Compact Remote)	(NC 2203, 2204)	£825
JBH2232		LiNX DLX-CR400LF (Compact Remote Low Force)	(NC 2203, 2204)	£1,075
JBH3000		mo-vis Micro LiNX	(NC 2201, 2203, 2204)	£2,960
JBH3001		mo-vis Multi LiNX	(NC 4401, 4402, 2201, 2203, 2204)	£2,640
JBH3110		mo-vis All Round Low Force Omni	(NC 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£2,200
JBH3111		mo-vis All Round Omni	(NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£1,600
JBC3112		mo-vis Heavy Duty Omni	(NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4300, 4301, 4302, 4400, 4401, 4402, 2201, 2203, 2204)	£3,000
COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)				
JBH2207		RH Mount	(NC 1707, 2209, 2213, 2217, 1404)	○
JBH2208		LH Mount	(NC 1708, 2210, 2214, 2218, 1405)	○
JBH2265		Centre Mount Swing Away to RH (Nucleus)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£495
JBH2266		Centre Mount Swing Away to LH (Nucleus)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£495
JBZ2215		Individual mounting option per request (Special). Requires additional information/specification.		POA
JBH3300		Maxx Resolve RH	(NC 4301, 4302)	£270
JBH3301		Maxx Resolve LH	(NC 4301, 4302)	£270
JBZ3310		BodyPoint Trilock Midline RH	(NC 4301, 4302, 0002)	£480
JBZ3311		BodyPoint Trilock Midline LH	(NC 4301, 4302, 0002)	£480
JBZ3320		Tray Mounted	(NC 4301, 4302, 0002)	○
JBZ3330		Individual mounting option per request (Special). Requires additional information/specification.		POA

CHIN CONTROLS				
JBH3000		mo-vis Micro LiNX	(NC 2201, 2203, 2204)	£2,960
JBH3001		mo-vis Multi LiNX	(NC 4401, 4402, 2201, 2203, 2204)	£2,640
JBH3110		mo-vis All Round Low Force Omni (consider option 3810 for switch access)	(NC 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£2,200
JBH3111		mo-vis All Round Omni (consider option 3810 for switch access)	(NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£1,600
CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)				
JBH3400		mo-vis Chin Control Harness (One Size)	(NC 4301, 4302)	£1,000
JBH3410		mo-vis Manual Lateral Swing Away RH Mounted	(NC 4301, 4302)	£1,520
JBH3411		mo-vis Manual Lateral Swing Away LH Mounted	(NC 4301, 4302)	£1,520
JBH3412		mo-vis Manual Vertical Swing Away RH Mounted	(NC 4301, 4302)	£1,520
JBH3413		mo-vis Manual Vertical Swing Away LH Mounted	(NC 4301, 4302)	£1,520
JBH3420		mo-vis Powered Lateral Swing Away RH Mounted	(NC 4301, 4302)	£2,400
JBH3421		mo-vis Powered Lateral Swing Away LH Mounted	(NC 4301, 4302)	£2,400
JBH3422		mo-vis Powered Vertical Swing Away RH Mounted	(NC 4301, 4302)	£2,400
JBH3423		mo-vis Powered Vertical Swing Away LH Mounted	(NC 4301, 4302)	£2,400
JBZ3330		Individual mounting option per request (Special). Requires additional information/specification.		POA
HEAD CONTROLS				
JBH3200		mo-vis Proportional Head Control Omni	(NC 1401, 2894, 2201, 2203, 2204)	£4,500
JBH3221		mo-vis Switched Head Control Omni	(NC 1401, 2894, 2201, 2203, 2204)	£4,500
JBH3222		mo-vis Switched Head Control w/ Sip & Puff Omni	(NC 1401, 2894, 2201, 2203, 2204)	£4,750
FOOT CONTROL				
JBH3010		mo-vis Foot Control LiNX RH Side	(NC 2203, 2204)	£3,300
JBH3011		mo-vis Foot Control LiNX LH Side	(NC 2203, 2204)	£3,300
LiNX 10-WAY SWITCH				
JBH2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access		£500
LiNX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)				
JBH2217		Bracket RH, mounted for user	(NC 2207, 2213, 2209, 2265, 1404)	○
JBH2218		Bracket LH, mounted for user	(NC 2208, 2214, 2210, 2266, 1405)	○
JBH2291		Bracket RH, mounted for attendant	(NC 2259, 1221, 2108)	£65
JBH2292		Bracket LH, mounted for attendant	(NC 2260, 1222, 2108)	£65
LiNX USB CHARGER PORT				
JBH2253		USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706)	£160
LiNX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)				
JBH2281		RH Mount		○
JBH2282		LH Mount		○
LiNX INPUT MODULE				
JBH2273		DLX-IN200 Module (Jack port for external power button)		£270
UNZ9249		DLX-IN500 Module (for DB9 Drive input. Only compatible with DLX-REM400 / DLX-REM500)	(NC 2201, 2203, 2203)	£785
FURTHER LiNX / MO-VIS COMPONENTS / ACCESSORIES				
JBH2302		Not Necessary	(NC 2254, 1711, 1712, 1713, 1714, 1715, 1716)	◆
JBH2300		To be enclosed only. Standard 300 mm Cable Length for mo-vis Switches		○
JBH2305		To be enclosed only. Extended 1500 mm Cable Length for mo-vis Switches		○
JBZ2301		To be connected and programmed (Special)		POA
BUDDY BUTTONS				
JBH1685		Buddy Button. Red / Yellow / Green / Blue / Black provided. Interchangeable colours.		£90

EGG SWITCH			
JBH2254		Egg Switch, black	£90
PIKO BUTTONS			
JBH1711		Red Piko Button	£160
JBH1712		Yellow Piko Button	£160
JBH1713		Green Piko Button	£160
JBH1714		Orange Piko Button	£160
JBH1715		Blue Piko Button	£160
JBH1716		Black Piko Button	£160
MO-VIS TWISTER BASIC			
JBH3500		Black	£140
JBH3501		Red	£140
JBH3502		Yellow	£140
JBH3503		Blue	£140
JBH3504		Green	£140
MO-VIS TWISTER BASIC MOUNTING (Only applicable for mo-vis Twister Basic)			
JBH3520		Mounting Set	£50
JBH3525		Table Mounting Set	£50
MO-VIS TWISTER			
JBH3600		Black	£160
JBH3601		Red	£160
JBH3602		Yellow	£160
JBH3603		Blue	£160
JBH3604		Green	£160
MO-VIS TWISTER ON BENDED TUBE INC. D6 CLAMP			
JBH3700		Black	£220
JBH3701		Red	£220
JBH3702		Yellow	£220
JBH3703		Blue	£220
JBH3704		Green	£220
MO-VIS SATELLITE TWISTER (Only applicable for mo-vis All Round Remotes)			
JBH3800		Black	£200
JBH3801		Red	£200
JBH3802		Yellow	£200
JBH3803		Blue	£200
JBH3804		Green	£200
JBH3810		All Round Chin Control Set (includes All Round Joystick Ball and 2x Satellite Switches)	£420
MO-VIS TWISTER GOOSENECK			
JBH3900		Black	£260
JBH3901		Red	£260
JBH3902		Green	£260

MO-VIS Q2M TWIST & SWING			
JBH4000		Black	£500
JBH4001		Red	£500
JBH4002		Yellow	£500
JBH4003		Blue	£500
JBH4004		Green	£500
MO-VIS TWISTER PRO			
JBH4100		Black	£160
JBH4101		Red	£160
JBH4102		Green	£160
MO-VIS MULTI SWITCH			
JBH4200		Multi Switch	£620
JBH4201		Multi Switch Proximity Sensors	£860
JBH4202		Multi Switch Proximity Set	£840
MO-VIS HANDWARMER			
JBH4300		Hand Warmer	£800
JBZ4301		Hand Warmer 2 Dome Hand Hood	£400
JBZ4302		Hand Warmer 1 Dome Hand Hood	£300
MO-VIS JOYSTICK KNOBS			
JBH4400		Micro Joystick Ball	£60
JBH4401		Micro Joystick Cup 10.1mm	£60
JBH4402		Micro Joystick Cup 15.5mm	£60
JBH4403		Multi Joystick Ball 4mm	£60
JBH4404		Multi Joystick Cup 24mm	£60
JBH4405		All Round Joystick Ball 6.35mm	£60
JBH4406		All Round Joystick Cup 37mm	£60
BODYPOINT JOYSTICK KNOBS			
JBH2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A	£75
JBH2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A	£75
JBH2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£75
JBH2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£75
5:1 SWITCH ADAPTOR			
JBH1681		5:1 Switch Adapter (9-pin sub D to 5 jack ports)	£120
SHROUD COLOUR*			
JBH2871		Mid range - Black/Anthracite Shroud	○
JBH2872		Mid range - Silver Shroud	○
JBH2873		Mid range - Bengal Red Shroud	○
JBH2874		Mid range - Mercury Blue Shroud	○
JBH2875		Mid range - Midnight Purple Shroud	○
JBZ0299		Premium - Individual Shroud Colour. RAL number required Extended lead lime: + 15 days.	£300

RIM COLOUR*			
JBH2882		Mid range - Black Rim Inserts	○
JBH2883		Mid range - Silver Rim Inserts	○
JBH2884		Mid range - Bengal Red Rim Inserts	○
JBH2885		Mid range - Mercury Blue Rim Inserts	○
JBH2886		Mid range - Midnight Purple Rim Inserts	○
JBZ2887		Premium - Individual Rim Insert Colour. RAL number required Extended lead lime: + 15 days. Price Included with Premium Shroud Kit	○
		(NC 2872, 2873, 2874, 2875)	
JBZ2887		Premium - Falcon White Rim Inserts	○
		(NC 2872, 2873, 2874, 2875, 0299)	
FORKS*			
JBH2631		Single sided fork	£250
JBH2632		Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆
CASTORS*			
JBH0835		9" (225 mm) x 2.80 (70 mm) Castor, pneumatic	○
JBH0836		9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	◆
CASTOR LOCKS			
JBH1490		Castor Lock	£195
DRIVE WHEELS*			
JBH0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
JBH0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	◆
C.T.C SUSPENSION*			
JBH2870		CTC Suspension Control - Traction - Comfort	◆
SUSPENSION SETTING*			
JBH1916		Suspension up to 160 kg User Weight For RX 40 HD only	◆
		(NC 0002)	
MOTOR TYPE*			
JBH1294		Indoor setting, reduced speed approx. 3 km/h	○
JBH1282		6 km/h	◆
JBH1908		10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights	£515
		(NC 2203, 1437, 1254)	
JBH1909		12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	£600
		(NC 2201, 2203, 1437, 1254)	
BATTERY TYPE*			
JBH1259		60 Ah Gel	◆
JBH1260		73.5 Ah Gel	£100
		(NC 1450)	
G TRAC			
JBH1215		G-Trac Module	£400
CHARGER*			
JBH1343		Opti Charge Charger 8A (IP44)	◆
JBH1313		Opti Charge Charger 10A (IP44)	£155
ANTI-TIPPERS*			
JBH1205		Anti-tipper	◆
LIGHTS			
JBH1437		LED Light system with indicators 30 Lux headlights, included with 10km/h motors	£450
		(NC 2203, 1908, 1909)	
JBH1254		Without Light system	○
		(NC 1909)	

KERB RIDER			
JBH1244		Kerb Rider with park position feature (NC 0556, 0589, 1294)	£215
TRANSPORTATION KIT*			
JBH1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.	£580
JBH1493		4-Point tie-down (For wheelchair transportation without user)	○
JBH1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. Not compatible with Seat lifter.	◆
JBH1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory. (NC 1091, 1092, 1093)	○
EASY ACCESS FRONT TIE DOWNS			
JBH2894		Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down. For transport worthy versions (1494/1401) only	£55
FIXATION BELT*			
JBH1329		Lap belt with Auto Lock Standard mount on backrest	○
JBH2113		Lap belt with Auto Lock Flexible mount in side rails	○
JBH2701		BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£95
JBH2702		BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm (NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£95
JBH2703		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£90
JBH2704		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm (NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£90
JBH2705		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£90
JBH2706		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm (NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£90
REAR VIEW MIRROR			
JBH1415		Rear view mirror RH (NC 2217, 2209, 2213, 2265)	£45
JBH1339		Rear view mirror LH (NC 2218, 2210, 2214, 2266)	£45
BAG HOOKS			
JBH2128		Bag Hooks (Max load up to 6kg = impact on max user weight) (NC 2172, 2173, 2227, 2228, 2291, 2292)	£80
CANE HOLDER			
JBH1312		Cane holder RH (NC 2227, 2259, 1221)	£45
JBH1363		Cane holder LH (NC 2228, 2260, 1222)	£45
WHEEL FENDERS			
JBH1416		Fender for Castors	£80
JBH1417		Fender for Drive wheels	£160
TRAY			
JBH1404		Tray, Swing Away RH (NC 1707, 2209, 2265, 2266, 2207, 2213, 2214, 2217)	£200
JBH1405		Tray, Swing Away LH (NC 1708, 2210, 2265, 2266, 2208, 2213, 2214, 2218)	£200
TOTAL PRICE			

Non Conforming Options:

JBH0002	AVIVA RX40 HD Modulite
JBH2112	Telescopic Seat Frame + Sling Seat with Side Rails. Width 380 mm to 530 mm. +30 mm seat height. Additional 50 mm width by arm rest adjustment.
JBH2121	380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.
JBH2122	430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.
JBH2123	480 mm Seat Plate
JBH2124	530 mm Seat Plate
JBH2104	550 mm Seat Plate
JBH2106	Body Fit Kit - Outrigging Backrest Brackets and Pivot Point (outrig: 1" /25mm per side)
JBH2132	410 mm
JBH2134	460 mm
JBH2142	Flo-tech Lite Cushion
JBH2147	Matrx Libra Cushion, 90 mm thick (Extended lead time)
JBH2144	Matrx Vi Cushion (Mesh fabric)
JBH2148	Standard Cushion (90mm high) - For 550mm Seat Plate Only
JBH2140	Without seat cushion
JBH1980	Tapered Wedges (Pair)
JBH1981	Abductor
JBH1982	Lateral Wedge
JBH1983	Anterior Seat Wedges
JBH1980	Tapered Wedges (Pair)
JBH1981	Abductor
JBH1982	Lateral Wedge
JBH1983	Anterior Seat Wedges
JBH1985	Pelvic Obliquity Kit (Contains 2x positioning accessories)
JBH2163	Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. Push Handles included.
JBH2108	Backposts without Push Handles
JBH2171	Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.
JBH2179	Flex3 Tension Adjustable Upholstery. Height 540 mm (seat plate to top of upholstery).
JBH2172	Flex3 adjustable Backrest Plates (without Backrest Cushion). Flex3 Backrest plate width and height adjustable 560 mm - 790 mm
JBH2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion. Flex3 Backrest plate width and height adjustable 560 mm - 790 mm
JBH2400	Support+ Standard Backrest (Spacer fabric)
JBH2401	Support+ Active Strapped Backrest (Spacer fabric)
JBH2402	Support+ Passive Strapped Backrest (Spacer fabric)
JBH2403	Support+ Extra Deep Backrest (Spacer fabric)
JBH2174	Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen
JBH2175	Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBH2176	Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBH2177	Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen

JBH2178	Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBH2168	Matrx MAC Back - Foam Cushion [Mesh fabric]
JBH2169	Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]
JBH2170	Without Backrest Upholstery (Back Posts/Back Frame only)
JBH4520	Tall - 570 mm (Measured from Seat Plate)
JBH4521	Short - 500 mm (Measured from Seat Plate)
JBH2380	405 mm (16 inch) MAC BH
JBH2381	455 mm (18 inch) MAC BH
JBH2382	505 mm (20 inch) MAC BH
JBH2383	560 mm (22 inch) MAC BH
JBH2385	90 mm (3.5 inch) MAC Contour Depth LH
JBH2386	150 mm (6 inch) MAC Contour Depth LH
JBH2387	230 mm (9 inch) MAC Contour Depth LH
JBH2389	90 mm (3.5 inch) MAC Contour Depth RH
JBH2390	150 mm (6 inch) MAC Contour Depth RH
JBH2391	230 mm (9 inch) MAC Contour Depth RH
JBH2393	Backrest -25 mm (-1 inch) to seat width
JBH2394	Backrest matches seat width
JBH2395	Backrest +25 mm (+1 inch) to seat width
JBH2151	Black Nylon Fabric TR18
JBH2152	Black Multi Stretch PU Coated Wipe Down Fabric TR26
JBH2155	Black Spacer Fabric TR36
JBH2156	Black Nylon (For 550mm SW Only)
JBH2181	T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm
JBH2198	Padded Armrest with padded side guard (Water Fall Type)
JBH2193	Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm).
JBH2192	Following Armrest. Height adjustable from 230 - 300 mm
JBH2194	Following Armrest. Height adjustable from 300 - 360 mm
JBH2195	Short Armpad (325 mm)
JBH2196	Long Armpad (375 mm)
JBH2197	Long Ergonomic Armpad (397 mm)
JBH2410	Support+ Standard Armpad Short (80 mm x 300 mm)
JBH2411	Support+ Standard Armpad Long (80 mm x 400 mm)
JBH2412	Support+ P-Shaped Armpad (120 mm x 400 mm)
JBH2199	Soft Armpad (REA Azalea Wide Pad)
JBH2184	Hip Support LH, pad size 140 x 100 mm
JBH2186	Hip Support LH, pad size 230 x 100 mm
JBH2188	Hip Support LH, pad size 330 x 100 mm
JBH2185	Hip Support RH, pad size 140 x 100 mm
JBH2187	Hip Support RH, pad size 230 x 100 mm
JBH2189	Hip Support RH, pad size 330 x 100 mm
JBH1350	Neckrest depth adjustable

JBH1432	Headrest sideways adjustable
JBH1427	Headrest with cheek support
JBH2328	Support+ Flex Headrest (Medium). Complete with Mounting Hardware & Stem (_2332).
JBH2329	Support+ Flex Headrest (Large). Complete with Mounting Hardware & Stem (_2332).
JBH2330	Support+ Standard Headrest (Medium). Complete with Mounting Hardware & Stem (_2332).
JBH2331	Support+ Standard Headrest (Large). Complete with Mounting Hardware & Stem (_2332).
JBH2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware
JBH2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware
JBH2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware
JBH1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
JBH1495	Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts
JBH4510	Support+ Wipedown Cushion Cover (PU Coated Type)
JBH4530	Support+ Wipedown Backrest Cover (PU Coated Type)
JBH4540	Support+ Wipedown Armpad Covers (PU Coated Type)
JBH4550	Support+ Wipedown Headrest Cover (PU Coated Type)
JBH1445	Lateral Support with extended arm, multi-adjustable pad
JBH1091	Manually settable Seat Tilt 0 - 15°
JBH1092	Powered "fixed pivot" Tilt 20°
JBH1093	Powered Tilt (0° - 30°) with CoG shift
JBH1094	Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift. 6 point Tie downs must be selected
JBH0977	Powered Recliner 30°
JBH0556	Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm
JBH0589	Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm
JBH0503	Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools
JBH0533	Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad
JBH0531	Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad
JBH0501	Without Leg Rests
JBH2365	Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")
JBH0559	Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")
JBH0516	Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5)
JBH0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
JBH0512	Fixed Footplate pair
JBH0504	Angle Adjustable Footplate pair
JBH0513	Angle and depth adjustable Footplate pair
JBH0501	Without footplates for option 0501 without Leg Rests
JBH2201	DLX-REM110 (Drive Only)
JBH2203	DLX-REM211 (Drive & Actuators)
JBH2204	DLX-REM216 (Drive, Actuators & Lights)
JBH2205	DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)
JBH2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
JBH2270	Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)

JBH1707	Swing Away RH Bracket
JBH1708	Swing Away LH Bracket
JBH1709	Matrx Resolve Swing in/out Bracket RH
JBH1710	Matrx Resolve Swing in/out Bracket LH
JBH2213	Centre Mount Swing Away to RH Nucleus Bracket
JBH2214	Centre Mount Swing Away to LH Nucleus Bracket
JBH2209	REM500 Swing Away RH Multi-Adjustable Bracket
JBH2210	REM500 Swing Away LH Multi-Adjustable Bracket
JBH2227	Bracket RH mounted for Attendant
JBH2228	Bracket LH mounted for Attendant
JBH2219	LiNX DLX-ACU200 (Compact Dual Control)
JBH2245	LiNX DLX-IDC (Intuitive Dual Control)
JBH2259	Bracket RH, Attendant
JBH2260	Bracket LH, Attendant
JBH1221	Bracket RH, detachable
JBH1222	Bracket LH, detachable
JBH9260	Bracket for mo-vis Scoot
JBH2231	LiNX DLX-CR400 (Compact Remote)
JBH2232	LiNX DLX-CR400LF (Compact Remote Low Force)
JBH2207	RH Mount
JBH2208	LH Mount
JBH2265	Centre Mount Swing Away to RH (Nucleus)
JBH2266	Centre Mount Swing Away to LH (Nucleus)
JBH3410	mo-vis Manual Lateral Swing Away RH Mounted
JBH3411	mo-vis Manual Lateral Swing Away LH Mounted
JBH3412	mo-vis Manual Vertical Swing Away RH Mounted
JBH3413	mo-vis Manual Vertical Swing Away LH Mounted
JBH3420	mo-vis Powered Lateral Swing Away RH Mounted
JBH3421	mo-vis Powered Lateral Swing Away LH Mounted
JBH3422	mo-vis Powered Vertical Swing Away RH Mounted
JBH3423	mo-vis Powered Vertical Swing Away LH Mounted
JBH3200	mo-vis Proportional Head Control Omni
JBH3221	mo-vis Switched Head Control Omni
JBH3222	mo-vis Switched Head Control w/ Sip & Puff Omni
JBH2217	Bracket RH, mounted for user
JBH2218	Bracket LH, mounted for user
JBH2291	Bracket RH, mounted for attendant
JBH2292	Bracket LH, mounted for attendant
JBH2253	USB Charger Port DLX-USB02-A
JBH2273	DLX-IN200 Module (Jack port for external power button)
JBH2254	Egg Switch, black
JBH1711	Red Piko Button
JBH1712	Yellow Piko Button
JBH1713	Green Piko Button

JBH1714	Orange Piko Button
JBH1715	Blue Piko Button
JBH1716	Black Piko Button
JBH4300	Hand Warmer
JBZ4301	Hand Warmer 2 Dome Hand Hood
JBZ4302	Hand Warmer 1 Dome Hand Hood
JBH4400	Micro Joystick Ball
JBH4401	Micro Joystick Cup 10.1mm
JBH4402	Micro Joystick Cup 15.5mm
JBH2872	Mid range - Silver Shroud
JBH2873	Mid range - Bengal Red Shroud
JBH2874	Mid range - Mercury Blue Shroud
JBH2875	Mid range - Midnight Purple Shroud
JBZ0299	Premium - Individual Shroud Colour. RAL number required Extended lead lime: + 15 days.
JBH1294	Indoor setting, reduced speed approx. 3 km/h
JBH1908	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBH1909	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBH1437	LED Light system with indicators 30 Lux headlights, included with 10km/h motors
JBH1254	Without Light system
JBH1244	Kerb Rider with park position feature
JBH1450	Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.
JBH1401	6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory.
JBH2894	Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down. For transport worthy versions (1494/1401) only
JBH2701	BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm
JBH2702	BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm
JBH2703	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm
JBH2704	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm
JBH2705	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm
JBH2706	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm
JBH1415	Rear view mirror RH
JBH1339	Rear view mirror LH
JBH2128	Bag Hooks (Max load up to 6kg = impact on max user weight)
JBH1312	Cane holder RH
JBH1363	Cane holder LH
JBH1416	Fender for Castors
JBH1404	Tray, Swing Away RH
JBH1405	Tray, Swing Away LH



For a more complex customer requirement,
please email: salesUK@dhcare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):