



**Yes, you can:**

Tel: 0800 043 0881  
Email: salesUK@dhcare.com  
Online Spares available at: www.dhcare.co.uk

# Power Retail Prescription Form

September 2026  
Customer Ref:

|                          |                    |
|--------------------------|--------------------|
| <b>Customer Ref*</b>     | <b>Account No*</b> |
| <b>Contact No*</b>       | <b>Order No*</b>   |
| <b>Delivery Address*</b> | <b>Order Qty*</b>  |
|                          | <b>Order date*</b> |

## AVIVA RX40 HD Modulite HD

### Standard Features

- Max user weight 160kg (25 Stone)
- Crash tested ISO7176 - 19
- Modulite HD Seating System
- HD telescopic adjustable seat widths from 530 to 550 mm with an additional 50mm by arms
- Seat depths from 510mm - 580mm, stepless adjustment
- Tension adjustable backrest
- Standard T-Armrests, height and depth adjustable, detachable, with short armpads
- Manual settable Seat Tilt (0° - 15°)
- Manual settable Backrest Recline (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplates
- LiNX DLX-REM110 with standard bracket
- Choice of 5 shroud colours and 5 rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension (Control - Traction - Comfort)
- 6 km/h (4mph) 4 pole motors and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £7,250**



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

**Code**                      **Description**

#### BASE MODEL\*

|         |                          |                                                                                                     |        |
|---------|--------------------------|-----------------------------------------------------------------------------------------------------|--------|
| JBH0002 | <input type="checkbox"/> | AVIVA RX40 HD Modulite HD                                                                           | £7,250 |
| JBZ1300 | <input type="checkbox"/> | Special AVIVA RX40 HD Modulite HD (Extended Lead Time). Please provide details in comments section. | POA    |

#### CHASSIS WIDTH\*

|         |                          |                           |   |
|---------|--------------------------|---------------------------|---|
| JBH2830 | <input type="checkbox"/> | Chassis Wide (635 mm/25") | ◆ |
|---------|--------------------------|---------------------------|---|

| SEATING TYPE*                                                        |  |                                                                                                                                                                                                |                                                                                                                                   |      |
|----------------------------------------------------------------------|--|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------|------|
| JBH2101                                                              |  | Modulite HD Seat Unit                                                                                                                                                                          | (NC 0002, 0002, 2273, 2149, 2400, 2401, 2402, 2403, 4520, 4521, 2410, 2411, 2412, 2328, 2329, 2330, 2331, 4510, 4530, 4540, 4550) | ○    |
| JBH2103                                                              |  | Modulite Support+                                                                                                                                                                              | (NC 2104, 2140, 2148, 2171, 2156, 2195, 2196, 2197, 1350, 1427, 1432, 1495, 2324, 2325, 2327)                                     | ○    |
| SEAT FRAME*                                                          |  |                                                                                                                                                                                                |                                                                                                                                   |      |
| JBH2102                                                              |  | HD Telescopic Seat Frame and Plate with Side Rails.<br>Width 530 to 550 mm                                                                                                                     |                                                                                                                                   | ○    |
| SEAT WIDTH*                                                          |  |                                                                                                                                                                                                |                                                                                                                                   |      |
| JBH2124                                                              |  | 530 mm Seat Plate                                                                                                                                                                              |                                                                                                                                   | ○    |
| JBH2104                                                              |  | 550 mm Seat Plate                                                                                                                                                                              | (NC 0556, 0589, 1404, 1405)                                                                                                       | ○    |
| SEAT DEPTH*                                                          |  |                                                                                                                                                                                                |                                                                                                                                   |      |
| JBH2136                                                              |  | 510 mm                                                                                                                                                                                         |                                                                                                                                   | ○    |
| JBH2137                                                              |  | 580 mm                                                                                                                                                                                         |                                                                                                                                   | ○    |
| SEAT CUSHION*                                                        |  |                                                                                                                                                                                                |                                                                                                                                   |      |
| JBH2148                                                              |  | Standard Cushion (90mm high) - For 550mm Seat Plate Only                                                                                                                                       |                                                                                                                                   | ○    |
| JBH2149                                                              |  | Support+ Cushion                                                                                                                                                                               |                                                                                                                                   | £270 |
| JBH2140                                                              |  | Without seat cushion                                                                                                                                                                           |                                                                                                                                   | ○    |
| BACKREST SUPPORT*                                                    |  |                                                                                                                                                                                                |                                                                                                                                   |      |
| JBH2162                                                              |  | Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/MatrX Backs, measured from seat plate to top of upholstery.<br>Width adjustable cast bracket. Push Handles included. |                                                                                                                                   | ◆    |
| BACKREST UPHOLSTERY / CUSHION*                                       |  |                                                                                                                                                                                                |                                                                                                                                   |      |
| JBH2171                                                              |  | Tension Adjustable Upholstery with Nylon overlay.<br>Short = 480 mm / 7 straps. Long = 540 mm / 8 straps.<br>Measured from seat plate to top of upholstery.                                    |                                                                                                                                   | ◆    |
| JBH2400                                                              |  | Support+ Standard Backrest (Spacer fabric)                                                                                                                                                     | (NC 2156)                                                                                                                         | £330 |
| JBH2401                                                              |  | Support+ Active Strapped Backrest (Spacer fabric)                                                                                                                                              | (NC 2156)                                                                                                                         | £400 |
| JBH2402                                                              |  | Support+ Passive Strapped Backrest (Spacer fabric)                                                                                                                                             | (NC 2156)                                                                                                                         | £400 |
| JBH2403                                                              |  | Support+ Extra Deep Backrest (Spacer fabric)                                                                                                                                                   | (NC 2156)                                                                                                                         | £450 |
| SUPPORT+ BACKREST HEIGHT (Only applicable for Support+ Back options) |  |                                                                                                                                                                                                |                                                                                                                                   |      |
| JBH4520                                                              |  | Tall - 570 mm (Measured from Seat Plate)                                                                                                                                                       |                                                                                                                                   | ○    |
| JBH4521                                                              |  | Short - 500 mm (Measured from Seat Plate)                                                                                                                                                      |                                                                                                                                   | ○    |
| UPHOLSTERY COLOUR*                                                   |  |                                                                                                                                                                                                |                                                                                                                                   |      |
| JBH2156                                                              |  | Black Nylon Fabric TR18                                                                                                                                                                        |                                                                                                                                   | ◆    |
| ARMRESTS*                                                            |  |                                                                                                                                                                                                |                                                                                                                                   |      |
| JBH2191                                                              |  | T-Arm detachable, Standard. Adjustment range from seat plate:<br>245 - 315 mm                                                                                                                  |                                                                                                                                   | ○    |
| JBH2181                                                              |  | T-Arm detachable, Extended height. Adjustment range from seat plate:<br>295 - 365 mm                                                                                                           |                                                                                                                                   | £40  |
| ARMPADS*                                                             |  |                                                                                                                                                                                                |                                                                                                                                   |      |
| JBH2195                                                              |  | Short Armpad (325 mm)                                                                                                                                                                          |                                                                                                                                   | ◆    |
| JBH2196                                                              |  | Long Armpad (375 mm)                                                                                                                                                                           |                                                                                                                                   | ○    |
| JBH2197                                                              |  | Long Ergonomic Armpad (397 mm)                                                                                                                                                                 |                                                                                                                                   | ○    |
| JBH2410                                                              |  | Support+ Standard Armpad Short (80 mm x 300 mm)                                                                                                                                                |                                                                                                                                   | £100 |
| JBH2411                                                              |  | Support+ Standard Armpad Long (80 mm x 400 mm)                                                                                                                                                 |                                                                                                                                   | £110 |
| JBH2412                                                              |  | Support+ P-Shaped Armpad (120 mm x 400 mm)                                                                                                                                                     |                                                                                                                                   | £130 |

| HEADREST / NECKREST (Adequate headrest option must be selected)* |  |                                                                                                                                    |        |
|------------------------------------------------------------------|--|------------------------------------------------------------------------------------------------------------------------------------|--------|
| JBH1350                                                          |  | Neckrest depth adjustable                                                                                                          | £80    |
| JBH1432                                                          |  | Headrest sideways adjustable                                                                                                       | £225   |
| JBH1427                                                          |  | Headrest with cheek support                                                                                                        | £225   |
| JBH2328                                                          |  | Support+ Flex Headrest (Medium). Complete with Mounting Hardware & Stem ( _ 2332).                                                 | £250   |
| JBH2329                                                          |  | Support+ Flex Headrest (Large). Complete with Mounting Hardware & Stem ( _ 2332).                                                  | £270   |
| JBH2330                                                          |  | Support+ Standard Headrest (Medium). Complete with Mounting Hardware & Stem ( _ 2332).                                             | £230   |
| JBH2331                                                          |  | Support+ Standard Headrest (Large). Complete with Mounting Hardware & Stem ( _ 2332).                                              | £250   |
| JBH2324                                                          |  | Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware                                 | £325   |
| JBH2325                                                          |  | Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware                               | £325   |
| JBH2327                                                          |  | Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware                              | £440   |
| JBH1495                                                          |  | Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts | ○      |
| SUPPORT+ ACCESSORIES                                             |  |                                                                                                                                    |        |
| JBH4510                                                          |  | Support+ Wipedown Cushion Cover (PU Coated Type)                                                                                   | £30    |
| JBH4530                                                          |  | Support+ Wipedown Backrest Cover (PU Coated Type)                                                                                  | £30    |
| JBH4540                                                          |  | Support+ Wipedown Armpad Covers (PU Coated Type)                                                                                   | £30    |
| JBH4550                                                          |  | Support+ Wipedown Headrest Cover (PU Coated Type)                                                                                  | £30    |
| MANUAL / POWERED SEATING OPTIONS*                                |  |                                                                                                                                    |        |
| JBH1091                                                          |  | Manually settable Seat Tilt 0 - 15° (NC 1401)                                                                                      | ○      |
| JBH1092                                                          |  | Powered "fixed pivot" Tilt 20° (NC 2201)                                                                                           | £400   |
| JBH2055                                                          |  | Powered "fixed pivot" Tilt 20° & Sit Smarter (NC 2201)                                                                             | £540   |
| JBH1093                                                          |  | Powered Tilt (0° - 30°) with CoG shift (NC 2201)                                                                                   | £1,100 |
| JBH2057                                                          |  | Powered Tilt (0° - 30°) with CoG shift & Sit Smarter (NC 2201)                                                                     | £1,490 |
| JBH1094                                                          |  | Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift. 6 point Tie downs must be selected (NC 2201, 2220)                     | £1,900 |
| JBH2058                                                          |  | Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift & Sit Smarter. 6 point Tie downs must be selected (NC 2201, 2220)       | £2,570 |
| SEAT TO GROUND HEIGHT*                                           |  |                                                                                                                                    |        |
| JBH2857                                                          |  | 435 mm                                                                                                                             | ○      |
| JBH2858                                                          |  | 460 mm                                                                                                                             | ○      |
| JBH2859                                                          |  | 485 mm                                                                                                                             | ○      |
| BACKREST RECLINER*                                               |  |                                                                                                                                    |        |
| JBH0975                                                          |  | Manual Settable Recliner (in 3.75° steps) 90° - 120°                                                                               | ◆      |
| JBH0977                                                          |  | Powered Recliner 30° (NC 2201)                                                                                                     | £500   |
| JBH0979                                                          |  | Powered Recliner 30° & Sit Smarter                                                                                                 | £670   |

| LEG REST TYPE*       |  |                                                                                                                          |                                                                         |
|----------------------|--|--------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------|
| JBH0556              |  | Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm                                       | (NC 0517, 0512, 0504, 0513, 1244, 1416, 0501) £370                      |
| JBH0589              |  | Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm                             | (NC 2201, 0512, 0504, 0513, 1244, 1416, 0501) £1,620                    |
| JBH0600              |  | Centre Mounted: LNX Powered Centre Mounted Legrest & Sit Smarter, length compensating. Length 340 - 410 mm               | (NC 2201, 0512, 0504, 0513, 1244, 1416, 0501) £2,190                    |
| JBH0503              |  | Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm.<br>Angle adjustment: 70 - 0° with tools                       | (NC 0559, 0516, 0517, 0501) ◆                                           |
| JBH0533              |  | Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm.<br>Angle adjustment: 70 - 0° with handwheel and calf pad      | (NC 2365, 0559, 0516, 0517, 0501) £235                                  |
| JBH0531              |  | Side Mounted: ADE pair. Lower leg length: 290 - 460 mm.<br>Powered angle adjustment: 80 - 0° with calf pad               | (NC 2201, 2365, 0559, 0516, 0517, 0501) £660                            |
| JBH0530              |  | Side Mounted: ADE pair & Sit Smarter. Lower leg length: 290 - 460 mm.<br>Powered angle adjustment: 80 - 0° with calf pad | (NC 2201, 2365, 0559, 0516, 0517, 0501) £900                            |
| JBH0501              |  | Without Leg Rests                                                                                                        | (NC 0504, 0513, 2365, 0559, 0516, 0517) ○                               |
| FOOTPLATE TYPE*      |  |                                                                                                                          |                                                                         |
| JBH2365              |  | Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")                                                       | £105                                                                    |
| JBH0559              |  | Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")                                                     | ○                                                                       |
| JBH0516              |  | Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5")                                                      | £135                                                                    |
| JBH0517              |  | Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")                                                         | £270                                                                    |
| JBH0512              |  | Fixed Footplate pair                                                                                                     | (NC 0533, 0531) ◆                                                       |
| JBH0504              |  | Angle Adjustable Footplate pair                                                                                          | £75                                                                     |
| JBH0513              |  | Angle and depth adjustable Footplate pair                                                                                | £120                                                                    |
| JBH0501              |  | Without footplates for option 0501 without Leg Rests                                                                     | (NC 0503, 0533, 0531) ○                                                 |
| LINX CONTROLS*       |  |                                                                                                                          |                                                                         |
| JBH2201              |  | DLX-REM110 (Drive Only)                                                                                                  | (NC 1092, 1093, 1094, 0977, 0589, 2231, 2232, 2219, 2245, 2270, 1908) ◆ |
| JBH2203              |  | DLX-REM211 (Drive & Actuators)                                                                                           | (NC 1437, 2231, 2232, 2270, 1908, 1909) ○                               |
| JBH2204              |  | DLX-REM216 (Drive, Actuators & Lights)                                                                                   | (NC 2232, 2270) £200                                                    |
| JBH2205              |  | DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)                                     | (NC 2220) £425                                                          |
| JBH2206              |  | DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).                           | (NC 2220, 2270, 1207, 1208, 1707, 1708, 1709, 1710) £765                |
| REM400 TOGGLE OPTION |  |                                                                                                                          |                                                                         |
| JBH2270              |  | Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)                                                                | (NC 2203, 2204, 2206) £315                                              |

| CONTROL BRACKETRY*                                                       |  |                                                                                                     |                                                                                                                             |        |
|--------------------------------------------------------------------------|--|-----------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|--------|
| JBH1207                                                                  |  | Standard Bracket RH                                                                                 | (NC 2207, 2265, 1404)                                                                                                       | ○      |
| JBH1208                                                                  |  | Standard Bracket LH                                                                                 | (NC 2208, 2266, 1405)                                                                                                       | ○      |
| JBH1707                                                                  |  | Swing Away RH Bracket                                                                               | (NC 2207, 2265, 1404)                                                                                                       | £135   |
| JBH1708                                                                  |  | Swing Away LH Bracket                                                                               | (NC 2208, 2266, 1405)                                                                                                       | £135   |
| JBH1709                                                                  |  | Matrx Resolve Swing in/out Bracket RH                                                               | (NC 2207, 2265, 1404)                                                                                                       | £245   |
| JBH1710                                                                  |  | Matrx Resolve Swing in/out Bracket LH                                                               | (NC 2208, 2266, 1405)                                                                                                       | £245   |
| JBH2213                                                                  |  | Centre Mount Swing Away to RH Nucleus Bracket                                                       | (NC 2207, 2265, 1415, 2217, 1404, 1405)                                                                                     | £495   |
| JBH2214                                                                  |  | Centre Mount Swing Away to LH Nucleus Bracket                                                       | (NC 2208, 2266, 1339, 2218, 1404, 1405)                                                                                     | £495   |
| JBH2209                                                                  |  | REM500 Swing Away RH Multi-Adjustable Bracket                                                       | (NC 2203, 2204, 2205, 2207, 2265, 1415, 2217, 1404)                                                                         | ○      |
| JBH2210                                                                  |  | REM500 Swing Away LH Multi-Adjustable Bracket                                                       | (NC 2203, 2204, 2205, 2208, 2266, 1339, 2218, 1405)                                                                         | ○      |
| JBH2227                                                                  |  | Bracket RH mounted for Attendant                                                                    | (NC 2220, 2128, 2219, 1221, 1312)                                                                                           | £65    |
| JBH2228                                                                  |  | Bracket LH mounted for Attendant                                                                    | (NC 2220, 2128, 2219, 1222, 1363)                                                                                           | £65    |
| JBZ2239                                                                  |  | Individual mounting option per request (Special).<br>Requires additional information/specification. |                                                                                                                             | POA    |
| DUAL CONTROLS                                                            |  |                                                                                                     |                                                                                                                             |        |
| JBH2220                                                                  |  | LiNX DLX-REM050 (Dual Control used as secondary remote)                                             | (NC 2273, 2205, 2206, 1221, 1222, 2227, 2228, 9260)                                                                         | £395   |
| JBH2219                                                                  |  | LiNX DLX-ACU200 (Compact Dual Control)                                                              | (NC 1221, 1222, 2227, 2228, 9260)                                                                                           | £395   |
| JBH2245                                                                  |  | LiNX DLX-IDC (Intuitive Dual Control)                                                               | (NC 2259, 2260, 9260)                                                                                                       | £745   |
| JBH2420                                                                  |  | mo-vis Scoot Control w/ dedicated bracket included.                                                 | (NC 2259, 2260, 2259, 2260, 1221, 1222)                                                                                     | £3,420 |
| DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)     |  |                                                                                                     |                                                                                                                             |        |
| JBH2259                                                                  |  | Dual Control Standard Bracket RH                                                                    | (NC 2291, 2245, 1312)                                                                                                       | ○      |
| JBH2260                                                                  |  | Dual Control Standard Bracket LH                                                                    | (NC 2292, 2245, 1363)                                                                                                       | ○      |
| JBH1221                                                                  |  | IDC Dual Control Detachable Bracket RH                                                              | (NC 2291, 2220, 2219, 2227, 1312)                                                                                           | ○      |
| JBH1222                                                                  |  | IDC Dual Control Detachable Bracket LH                                                              | (NC 2292, 2220, 2219, 2228, 1363)                                                                                           | ○      |
| JBH9260                                                                  |  | Bracket for mo-vis Scoot                                                                            | (NC 2219, 2245)                                                                                                             | ○      |
| COMPACT REMOTES                                                          |  |                                                                                                     |                                                                                                                             |        |
| JBH2231                                                                  |  | LiNX DLX-CR400 (Compact Remote)                                                                     | (NC 2203, 2204)                                                                                                             | £825   |
| JBH2232                                                                  |  | LiNX DLX-CR400LF (Compact Remote Low Force)                                                         | (NC 2203, 2204)                                                                                                             | £1,075 |
| JBH3000                                                                  |  | mo-vis Micro LiNX                                                                                   | (NC 2201, 2203, 2204)                                                                                                       | £2,960 |
| JBH3001                                                                  |  | mo-vis Multi LiNX                                                                                   | (NC 4401, 4402, 2201, 2203, 2204)                                                                                           | £2,640 |
| JBH3110                                                                  |  | mo-vis All Round Low Force Omni                                                                     | (NC 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)                                                                         | £2,200 |
| JBH3111                                                                  |  | mo-vis All Round Omni                                                                               | (NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)                   | £1,600 |
| JBH3112                                                                  |  | mo-vis Heavy Duty Omni                                                                              | (NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4300, 4301, 4302, 4400, 4401, 4402, 2201, 2203, 2204) | £3,000 |
| COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected) |  |                                                                                                     |                                                                                                                             |        |
| JBH2207                                                                  |  | RH Mount                                                                                            | (NC 1707, 2209, 2213, 2217, 1404)                                                                                           | ○      |
| JBH2208                                                                  |  | LH Mount                                                                                            | (NC 1708, 2210, 2214, 2218, 1405)                                                                                           | ○      |
| JBH2265                                                                  |  | Centre Mount Swing Away to RH (Nucleus)                                                             | (NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)                                                                               | £495   |
| JBH2266                                                                  |  | Centre Mount Swing Away to LH (Nucleus)                                                             | (NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)                                                                               | £495   |
| JBZ2215                                                                  |  | Individual mounting option per request (Special).<br>Requires additional information/specification. |                                                                                                                             | POA    |
| JBH3300                                                                  |  | Maxx Resolve RH                                                                                     | (NC 4301, 4302)                                                                                                             | £270   |
| JBH3301                                                                  |  | Maxx Resolve LH                                                                                     | (NC 4301, 4302)                                                                                                             | £270   |
| JBZ3310                                                                  |  | BodyPoint Trilock Midline RH                                                                        | (NC 4301, 4302, 0002)                                                                                                       | £480   |
| JBZ3311                                                                  |  | BodyPoint Trilock Midline LH                                                                        | (NC 4301, 4302, 0002)                                                                                                       | £480   |
| JBZ3320                                                                  |  | Tray Mounted                                                                                        | (NC 4301, 4302, 0002)                                                                                                       | ○      |
| JBZ3330                                                                  |  | Individual mounting option per request (Special).<br>Requires additional information/specification. |                                                                                                                             | POA    |

| CHIN CONTROLS                                                                |  |                                                                                                  |                                                                                                           |        |
|------------------------------------------------------------------------------|--|--------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------|--------|
| JBH3000                                                                      |  | mo-vis Micro LiNX                                                                                | (NC 2201, 2203, 2204)                                                                                     | £2,960 |
| JBH3001                                                                      |  | mo-vis Multi LiNX                                                                                | (NC 4401, 4402, 2201, 2203, 2204)                                                                         | £2,640 |
| JBH3110                                                                      |  | mo-vis All Round Low Force Omni (consider option 3810 for switch access)                         | (NC 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)                                                       | £2,200 |
| JBH3111                                                                      |  | mo-vis All Round Omni (consider option 3810 for switch access)                                   | (NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204) | £1,600 |
| CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)         |  |                                                                                                  |                                                                                                           |        |
| JBH3400                                                                      |  | mo-vis Chin Control Harness (One Size)                                                           | (NC 4301, 4302)                                                                                           | £1,000 |
| JBH3410                                                                      |  | mo-vis Manual Lateral Swing Away RH Mounted                                                      | (NC 4301, 4302)                                                                                           | £1,520 |
| JBH3411                                                                      |  | mo-vis Manual Lateral Swing Away LH Mounted                                                      | (NC 4301, 4302)                                                                                           | £1,520 |
| JBH3412                                                                      |  | mo-vis Manual Vertical Swing Away RH Mounted                                                     | (NC 4301, 4302)                                                                                           | £1,520 |
| JBH3413                                                                      |  | mo-vis Manual Vertical Swing Away LH Mounted                                                     | (NC 4301, 4302)                                                                                           | £1,520 |
| JBH3420                                                                      |  | mo-vis Powered Lateral Swing Away RH Mounted                                                     | (NC 4301, 4302)                                                                                           | £2,400 |
| JBH3421                                                                      |  | mo-vis Powered Lateral Swing Away LH Mounted                                                     | (NC 4301, 4302)                                                                                           | £2,400 |
| JBH3422                                                                      |  | mo-vis Powered Vertical Swing Away RH Mounted                                                    | (NC 4301, 4302)                                                                                           | £2,400 |
| JBH3423                                                                      |  | mo-vis Powered Vertical Swing Away LH Mounted                                                    | (NC 4301, 4302)                                                                                           | £2,400 |
| JBZ3330                                                                      |  | Individual mounting option per request (Special). Requires additional information/specification. |                                                                                                           | POA    |
| HEAD CONTROLS                                                                |  |                                                                                                  |                                                                                                           |        |
| JBH3200                                                                      |  | mo-vis Proportional Head Control Omni                                                            | (NC 1401, 2894, 2201, 2203, 2204)                                                                         | £4,500 |
| JBH3221                                                                      |  | mo-vis Switched Head Control Omni                                                                | (NC 1401, 2894, 2201, 2203, 2204)                                                                         | £4,500 |
| JBH3222                                                                      |  | mo-vis Switched Head Control w/ Sip & Puff Omni                                                  | (NC 1401, 2894, 2201, 2203, 2204)                                                                         | £4,750 |
| FOOT CONTROL                                                                 |  |                                                                                                  |                                                                                                           |        |
| JBH3010                                                                      |  | mo-vis Foot Control LiNX RH Side                                                                 | (NC 2203, 2204, 2220)                                                                                     | £3,300 |
| JBH3011                                                                      |  | mo-vis Foot Control LiNX LH Side                                                                 | (NC 2203, 2204, 2220)                                                                                     | £3,300 |
| LiNX 10-WAY SWITCH                                                           |  |                                                                                                  |                                                                                                           |        |
| JBH2251                                                                      |  | 10-Way Switch DLX-SW10-A for direct Powered Seating Access                                       |                                                                                                           | £500   |
| LiNX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)  |  |                                                                                                  |                                                                                                           |        |
| JBH2217                                                                      |  | Bracket RH, mounted for user                                                                     | (NC 2207, 2213, 2209, 2265, 1404)                                                                         | ○      |
| JBH2218                                                                      |  | Bracket LH, mounted for user                                                                     | (NC 2208, 2214, 2210, 2266, 1405)                                                                         | ○      |
| JBH2291                                                                      |  | Bracket RH, mounted for attendant                                                                | (NC 2259, 1221)                                                                                           | £65    |
| JBH2292                                                                      |  | Bracket LH, mounted for attendant                                                                | (NC 2260, 1222)                                                                                           | £65    |
| LiNX USB CHARGER PORT                                                        |  |                                                                                                  |                                                                                                           |        |
| JBH2253                                                                      |  | USB Charger Port DLX-USB02-A                                                                     |                                                                                                           | £160   |
| LiNX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected) |  |                                                                                                  |                                                                                                           |        |
| JBH2281                                                                      |  | RH Mount                                                                                         |                                                                                                           | ○      |
| JBH2282                                                                      |  | LH Mount                                                                                         |                                                                                                           | ○      |
| LiNX INPUT MODULE                                                            |  |                                                                                                  |                                                                                                           |        |
| JBH2273                                                                      |  | DLX-IN200 Module (Jack port for external power button)                                           |                                                                                                           | £270   |
| UNZ9249                                                                      |  | DLX-IN500 Module (for DB9 Drive input. Only compatible with DLX-REM400 / DLX-REM500)             | (NC 2201, 2203, 2203)                                                                                     | £785   |
| FURTHER LiNX / MO-VIS COMPONENTS / ACCESSORIES                               |  |                                                                                                  |                                                                                                           |        |
| JBH2302                                                                      |  | Not Necessary                                                                                    | (NC 2254, 1711, 1712, 1713, 1714, 1715, 1716)                                                             | ◆      |
| JBH2300                                                                      |  | To be enclosed only. Standard 300 mm Cable Length for mo-vis Switches                            |                                                                                                           | ○      |
| JBH2305                                                                      |  | To be enclosed only. Extended 1500 mm Cable Length for mo-vis Switches                           |                                                                                                           | ○      |
| JBZ2301                                                                      |  | To be connected and programmed (Special)                                                         |                                                                                                           | POA    |
| BUDDY BUTTONS                                                                |  |                                                                                                  |                                                                                                           |        |
| JBH1685                                                                      |  | Buddy Button. Red / Yellow / Green / Blue / Black provided. Interchangeable colours.             |                                                                                                           | £90    |

| EGG SWITCH                                                               |  |                                                                                         |      |
|--------------------------------------------------------------------------|--|-----------------------------------------------------------------------------------------|------|
| JBH2254                                                                  |  | Egg Switch, black                                                                       | £90  |
| PIKO BUTTONS                                                             |  |                                                                                         |      |
| JBH1711                                                                  |  | Red Piko Button                                                                         | £160 |
| JBH1712                                                                  |  | Yellow Piko Button                                                                      | £160 |
| JBH1713                                                                  |  | Green Piko Button                                                                       | £160 |
| JBH1714                                                                  |  | Orange Piko Button                                                                      | £160 |
| JBH1715                                                                  |  | Blue Piko Button                                                                        | £160 |
| JBH1716                                                                  |  | Black Piko Button                                                                       | £160 |
| MO-VIS TWISTER BASIC                                                     |  |                                                                                         |      |
| JBH3500                                                                  |  | Black                                                                                   | £140 |
| JBH3501                                                                  |  | Red                                                                                     | £140 |
| JBH3502                                                                  |  | Yellow                                                                                  | £140 |
| JBH3503                                                                  |  | Blue                                                                                    | £140 |
| JBH3504                                                                  |  | Green                                                                                   | £140 |
| MO-VIS TWISTER BASIC MOUNTING (Only applicable for mo-vis Twister Basic) |  |                                                                                         |      |
| JBH3520                                                                  |  | Mounting Set                                                                            | £50  |
| JBH3525                                                                  |  | Table Mounting Set                                                                      | £50  |
| MO-VIS TWISTER                                                           |  |                                                                                         |      |
| JBH3600                                                                  |  | Black                                                                                   | £160 |
| JBH3601                                                                  |  | Red                                                                                     | £160 |
| JBH3602                                                                  |  | Yellow                                                                                  | £160 |
| JBH3603                                                                  |  | Blue                                                                                    | £160 |
| JBH3604                                                                  |  | Green                                                                                   | £160 |
| MO-VIS TWISTER ON BENDED TUBE INC. D6 CLAMP                              |  |                                                                                         |      |
| JBH3700                                                                  |  | Black                                                                                   | £220 |
| JBH3701                                                                  |  | Red                                                                                     | £220 |
| JBH3702                                                                  |  | Yellow                                                                                  | £220 |
| JBH3703                                                                  |  | Blue                                                                                    | £220 |
| JBH3704                                                                  |  | Green                                                                                   | £220 |
| MO-VIS SATELLITE TWISTER (Only applicable for mo-vis All Round Remotes)  |  |                                                                                         |      |
| JBH3800                                                                  |  | Black                                                                                   | £200 |
| JBH3801                                                                  |  | Red                                                                                     | £200 |
| JBH3802                                                                  |  | Yellow                                                                                  | £200 |
| JBH3803                                                                  |  | Blue                                                                                    | £200 |
| JBH3804                                                                  |  | Green                                                                                   | £200 |
| JBH3810                                                                  |  | All Round Chin Control Set (includes All Round Joystick Ball and 2x Satellite Switches) | £420 |
| MO-VIS TWISTER GOOSENECK                                                 |  |                                                                                         |      |
| JBH3900                                                                  |  | Black                                                                                   | £260 |
| JBH3901                                                                  |  | Red                                                                                     | £260 |
| JBH3902                                                                  |  | Green                                                                                   | £260 |

| MO-VIS Q2M TWIST & SWING |  |                                                                                           |      |
|--------------------------|--|-------------------------------------------------------------------------------------------|------|
| JBH4000                  |  | Black                                                                                     | £500 |
| JBH4001                  |  | Red                                                                                       | £500 |
| JBH4002                  |  | Yellow                                                                                    | £500 |
| JBH4003                  |  | Blue                                                                                      | £500 |
| JBH4004                  |  | Green                                                                                     | £500 |
| MO-VIS TWISTER PRO       |  |                                                                                           |      |
| JBH4100                  |  | Black                                                                                     | £160 |
| JBH4101                  |  | Red                                                                                       | £160 |
| JBH4102                  |  | Green                                                                                     | £160 |
| MO-VIS MULTI SWITCH      |  |                                                                                           |      |
| JBH4200                  |  | Multi Switch                                                                              | £620 |
| JBH4201                  |  | Multi Switch Proximity Sensors                                                            | £860 |
| JBH4202                  |  | Multi Switch Proximity Set                                                                | £840 |
| MO-VIS HANDWARMER        |  |                                                                                           |      |
| JBH4300                  |  | Hand Warmer                                                                               | £800 |
| JBZ4301                  |  | Hand Warmer 2 Dome Hand Hood                                                              | £400 |
| JBZ4302                  |  | Hand Warmer 1 Dome Hand Hood                                                              | £300 |
| MO-VIS JOYSTICK KNOBS    |  |                                                                                           |      |
| JBH4400                  |  | Micro Joystick Ball                                                                       | £60  |
| JBH4401                  |  | Micro Joystick Cup 10.1mm                                                                 | £60  |
| JBH4402                  |  | Micro Joystick Cup 15.5mm                                                                 | £60  |
| JBH4403                  |  | Multi Joystick Ball 4mm                                                                   | £60  |
| JBH4404                  |  | Multi Joystick Cup 24mm                                                                   | £60  |
| JBH4405                  |  | All Round Joystick Ball 6.35mm                                                            | £60  |
| JBH4406                  |  | All Round Joystick Cup 37mm                                                               | £60  |
| BODYPOINT JOYSTICK KNOBS |  |                                                                                           |      |
| JBH2721                  |  | 3 inch U-Shaped, Joystick Handle, Bodypoint PC101A                                        | £75  |
| JBH2722                  |  | 4 inch U-Shaped, Joystick Handle, Bodypoint PC102A                                        | £75  |
| JBH2723                  |  | Domed Rubber, Joystick Handle, Bodypoint PC107A                                           | £75  |
| JBH2724                  |  | Mushroom Grooved, Joystick Handle, Bodypoint PC110A                                       | £75  |
| 5:1 SWITCH ADAPTOR       |  |                                                                                           |      |
| JBH1681                  |  | 5:1 Switch Adapter ( 9-pin sub D to 5 jack ports)                                         | £120 |
| SHROUD COLOUR*           |  |                                                                                           |      |
| JBH2871                  |  | Mid range - Black/Anthracite Shroud                                                       | ○    |
| JBH2872                  |  | Mid range - Silver Shroud                                                                 | ○    |
| JBH2873                  |  | Mid range - Bengal Red Shroud                                                             | ○    |
| JBH2874                  |  | Mid range - Mercury Blue Shroud                                                           | ○    |
| JBH2875                  |  | Mid range - Midnight Purple Shroud                                                        | ○    |
| JBZ0299                  |  | Premium - Individual Shroud Colour. RAL number required<br>Extended lead lime: + 15 days. | £300 |

| RIM COLOUR*         |  |                                                                                                                                          |      |
|---------------------|--|------------------------------------------------------------------------------------------------------------------------------------------|------|
| JBH2882             |  | Mid range - Black Rim Inserts                                                                                                            | ○    |
| JBH2883             |  | Mid range - Silver Rim Inserts                                                                                                           | ○    |
| JBH2884             |  | Mid range - Bengal Red Rim Inserts                                                                                                       | ○    |
| JBH2885             |  | Mid range - Mercury Blue Rim Inserts                                                                                                     | ○    |
| JBH2886             |  | Mid range - Midnight Purple Rim Inserts                                                                                                  | ○    |
| JBZ2887             |  | Premium - Individual Rim Insert Colour. RAL number required<br>Extended lead lime: + 15 days. Price Included with Premium Shroud Kit     | ○    |
| FORKS*              |  |                                                                                                                                          |      |
| JBH2631             |  | Single sided fork                                                                                                                        | £250 |
| JBH2632             |  | Double sided fork Without castor rim insert<br>(enlarges wide base to 656 mm)                                                            | ◆    |
| CASTORS*            |  |                                                                                                                                          |      |
| JBH0835             |  | 9" (225 mm) x 2.80 (70 mm) Castor, pneumatic                                                                                             | ○    |
| JBH0836             |  | 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre                                                                             | ◆    |
| CASTOR LOCKS        |  |                                                                                                                                          |      |
| JBH1490             |  | Castor Lock                                                                                                                              | £195 |
| DRIVE WHEELS*       |  |                                                                                                                                          |      |
| JBH0704             |  | DW 14" (355 mm) x 3" (75 mm), Pneumatic                                                                                                  | ○    |
| JBH0705             |  | DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert                                                                       | ◆    |
| C.T.C SUSPENSION*   |  |                                                                                                                                          |      |
| JBH2870             |  | CTC Suspension Control - Traction - Comfort                                                                                              | ◆    |
| SUSPENSION SETTING* |  |                                                                                                                                          |      |
| JBH1916             |  | Suspension up to 160 kg User Weight For RX 40 HD (NC 0002) only                                                                          | ◆    |
| MOTOR TYPE*         |  |                                                                                                                                          |      |
| JBH1294             |  | Indoor setting, reduced speed approx. 3 km/h                                                                                             | ○    |
| JBH1282             |  | 6 km/h                                                                                                                                   | ◆    |
| JBH1908             |  | 10 km/h Package<br>Always combined with LED Light system with Indicators, 30 Lux headlights (NC 2203, 1437, 1254)                        | £515 |
| JBH1909             |  | 12 km/h Ultra High Speed Package<br>Always combined with LED Light system with Indicators, 30 Lux headlights (NC 2201, 2203, 1437, 1254) | £600 |
| BATTERY TYPE*       |  |                                                                                                                                          |      |
| JBH1259             |  | 60 Ah Gel                                                                                                                                | ◆    |
| JBH1260             |  | 73.5 Ah Gel (NC 1450)                                                                                                                    | £100 |
| G TRAC              |  |                                                                                                                                          |      |
| JBH1215             |  | G-Trac Module                                                                                                                            | £400 |
| CHARGER*            |  |                                                                                                                                          |      |
| JBH1343             |  | Opti Charge Charger 8A (IP44)                                                                                                            | ◆    |
| JBH1313             |  | Opti Charge Charger 10A (IP44)                                                                                                           | £155 |
| ANTI-TIPPERS*       |  |                                                                                                                                          |      |
| JBH1205             |  | Anti-tipper                                                                                                                              | ◆    |
| LIGHTS              |  |                                                                                                                                          |      |
| JBH1437             |  | LED Light system with indicators 30 Lux headlights, included with 10km/h motors (NC 2203, 1908, 1909)                                    | £450 |
| JBH1254             |  | Without Light system (NC 1909)                                                                                                           | ○    |

| KERB RIDER                  |  |                                                                                                                                                                                                                    |      |
|-----------------------------|--|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------|
| JBH1244                     |  | Kerb Rider with park position feature (NC 0556, 0589, 1294)                                                                                                                                                        | £215 |
| TRANSPORTATION KIT*         |  |                                                                                                                                                                                                                    |      |
| JBH1450                     |  | Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering. | £580 |
| JBH1493                     |  | 4-Point tie-down (For wheelchair transportation without user)                                                                                                                                                      | ○    |
| JBH1494                     |  | 4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. Not compatible with Seat lifter.                                               | ◆    |
| JBH1401                     |  | 6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory. (NC 1091, 1092, 1093)  | ○    |
| EASY ACCESS FRONT TIE DOWNS |  |                                                                                                                                                                                                                    |      |
| JBH2894                     |  | Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down. For transport worthy versions (1494/1401) only                                                                          | £55  |
| FIXATION BELT*              |  |                                                                                                                                                                                                                    |      |
| JBH1329                     |  | Lap belt with Auto Lock Standard mount on backrest                                                                                                                                                                 | ○    |
| JBH2113                     |  | Lap belt with Auto Lock Flexible mount in side rails                                                                                                                                                               | ○    |
| REAR VIEW MIRROR            |  |                                                                                                                                                                                                                    |      |
| JBH1415                     |  | Rear view mirror RH (NC 2217, 2209, 2213, 2265)                                                                                                                                                                    | £45  |
| JBH1339                     |  | Rear view mirror LH (NC 2218, 2210, 2214, 2266)                                                                                                                                                                    | £45  |
| BAG HOOKS                   |  |                                                                                                                                                                                                                    |      |
| JBH2128                     |  | Bag Hooks (Max load up to 6kg = impact on max user weight) (NC 2227, 2228, 2291, 2292)                                                                                                                             | £80  |
| CANE HOLDER                 |  |                                                                                                                                                                                                                    |      |
| JBH1312                     |  | Cane holder RH (NC 2227, 2259, 1221)                                                                                                                                                                               | £45  |
| JBH1363                     |  | Cane holder LH (NC 2228, 2260, 1222)                                                                                                                                                                               | £45  |
| WHEEL FENDERS               |  |                                                                                                                                                                                                                    |      |
| JBH1416                     |  | Fender for Castors                                                                                                                                                                                                 | £80  |
| JBH1417                     |  | Fender for Drive wheels                                                                                                                                                                                            | £160 |
| TRAY                        |  |                                                                                                                                                                                                                    |      |
| JBH1404                     |  | Tray, Swing Away RH (NC 1707, 2209, 2265, 2266, 2207, 2213, 2214, 2217)                                                                                                                                            | £200 |
| JBH1405                     |  | Tray, Swing Away LH (NC 1708, 2210, 2265, 2266, 2208, 2213, 2214, 2218)                                                                                                                                            | £200 |

|             |  |  |  |
|-------------|--|--|--|
| TOTAL PRICE |  |  |  |
|-------------|--|--|--|

| Non Conforming Options: |                                                                                                                                                             |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| JBH0002                 | AVIVA RX40 HD Modulite HD                                                                                                                                   |
| JBH2104                 | 550 mm Seat Plate                                                                                                                                           |
| JBH2148                 | Standard Cushion (90mm high) - For 550mm Seat Plate Only                                                                                                    |
| JBH2149                 | Support+ Cushion                                                                                                                                            |
| JBH2140                 | Without seat cushion                                                                                                                                        |
| JBH2171                 | Tension Adjustable Upholstery with Nylon overlay.<br>Short = 480 mm / 7 straps. Long = 540 mm / 8 straps.<br>Measured from seat plate to top of upholstery. |
| JBH2400                 | Support+ Standard Backrest (Spacer fabric)                                                                                                                  |
| JBH2401                 | Support+ Active Strapped Backrest (Spacer fabric)                                                                                                           |
| JBH2402                 | Support+ Passive Strapped Backrest (Spacer fabric)                                                                                                          |
| JBH2403                 | Support+ Extra Deep Backrest (Spacer fabric)                                                                                                                |
| JBH4520                 | Tall - 570 mm (Measured from Seat Plate)                                                                                                                    |
| JBH4521                 | Short - 500 mm (Measured from Seat Plate)                                                                                                                   |
| JBH2156                 | Black Nylon Fabric TR18                                                                                                                                     |
| JBH2195                 | Short Armpad (325 mm)                                                                                                                                       |
| JBH2196                 | Long Armpad (375 mm)                                                                                                                                        |
| JBH2197                 | Long Ergonomic Armpad (397 mm)                                                                                                                              |
| JBH2410                 | Support+ Standard Armpad Short (80 mm x 300 mm)                                                                                                             |
| JBH2411                 | Support+ Standard Armpad Long (80 mm x 400 mm)                                                                                                              |
| JBH2412                 | Support+ P-Shaped Armpad (120 mm x 400 mm)                                                                                                                  |
| JBH1350                 | Neckrest depth adjustable                                                                                                                                   |
| JBH1432                 | Headrest sideways adjustable                                                                                                                                |
| JBH1427                 | Headrest with cheek support                                                                                                                                 |
| JBH2328                 | Support+ Flex Headrest (Medium). Complete with Mounting Hardware & Stem ( _2332).                                                                           |
| JBH2329                 | Support+ Flex Headrest (Large). Complete with Mounting Hardware & Stem ( _2332).                                                                            |
| JBH2330                 | Support+ Standard Headrest (Medium). Complete with Mounting Hardware & Stem ( _2332).                                                                       |
| JBH2331                 | Support+ Standard Headrest (Large). Complete with Mounting Hardware & Stem ( _2332).                                                                        |
| JBH2324                 | Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware                                                          |
| JBH2325                 | Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware                                                        |
| JBH2327                 | Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware                                                       |
| JBH1495                 | Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts                          |
| JBH4510                 | Support+ Wipedown Cushion Cover (PU Coated Type)                                                                                                            |
| JBH4530                 | Support+ Wipedown Backrest Cover (PU Coated Type)                                                                                                           |
| JBH4540                 | Support+ Wipedown Armpad Covers (PU Coated Type)                                                                                                            |
| JBH4550                 | Support+ Wipedown Headrest Cover (PU Coated Type)                                                                                                           |
| JBH1091                 | Manually settable Seat Tilt 0 - 15°                                                                                                                         |
| JBH1092                 | Powered "fixed pivot" Tilt 20°                                                                                                                              |
| JBH1093                 | Powered Tilt (0° - 30°) with CoG shift                                                                                                                      |
| JBH1094                 | Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift. 6 point Tie downs must be selected                                                              |
| JBH0977                 | Powered Recliner 30°                                                                                                                                        |
| JBH0556                 | Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm                                                                          |
| JBH0589                 | Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm                                                                |

|         |                                                                                                                     |
|---------|---------------------------------------------------------------------------------------------------------------------|
| JBH0503 | Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm.<br>Angle adjustment: 70 - 0° with tools                  |
| JBH0533 | Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm.<br>Angle adjustment: 70 - 0° with handwheel and calf pad |
| JBH0531 | Side Mounted: ADE pair. Lower leg length: 290 - 460 mm.<br>Powered angle adjustment: 80 - 0° with calf pad          |
| JBH0501 | Without Leg Rests                                                                                                   |
| JBH2365 | Angle & width adjustable footplates S<br>(127 W x 190 mm D / 5"x7.5")                                               |
| JBH0559 | Angle & width adjustable footplates M<br>(140 W x 240 mm D / 5.5"x9.5")                                             |
| JBH0516 | Angle & width adjustable footplates L<br>(150 W x 290 mm D / 6"x11.5")                                              |
| JBH0517 | Angle adjustable foot platform M<br>(305 W x 275 mm D / 12"x10.75")                                                 |
| JBH0512 | Fixed Footplate pair                                                                                                |
| JBH0504 | Angle Adjustable Footplate pair                                                                                     |
| JBH0513 | Angle and depth adjustable Footplate pair                                                                           |
| JBH0501 | Without footplates for option 0501 without Leg Rests                                                                |
| JBH2201 | DLX-REM110 (Drive Only)                                                                                             |
| JBH2203 | DLX-REM211 (Drive & Actuators)                                                                                      |
| JBH2204 | DLX-REM216 (Drive, Actuators & Lights)                                                                              |
| JBH2205 | DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)                                |
| JBH2206 | DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).                      |
| JBH2270 | Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)                                                           |
| JBH1207 | Standard Bracket RH                                                                                                 |
| JBH1208 | Standard Bracket LH                                                                                                 |
| JBH1707 | Swing Away RH Bracket                                                                                               |
| JBH1708 | Swing Away LH Bracket                                                                                               |
| JBH1709 | Matrx Resolve Swing in/out Bracket RH                                                                               |
| JBH1710 | Matrx Resolve Swing in/out Bracket LH                                                                               |
| JBH2213 | Centre Mount Swing Away to RH Nucleus Bracket                                                                       |
| JBH2214 | Centre Mount Swing Away to LH Nucleus Bracket                                                                       |
| JBH2209 | REM500 Swing Away RH Multi-Adjustable Bracket                                                                       |
| JBH2210 | REM500 Swing Away LH Multi-Adjustable Bracket                                                                       |
| JBH2227 | Bracket RH mounted for Attendant                                                                                    |
| JBH2228 | Bracket LH mounted for Attendant                                                                                    |
| JBH2220 | LiNX DLX-REM050 (Dual Control used as secondary remote)                                                             |
| JBH2219 | LiNX DLX-ACU200 (Compact Dual Control)                                                                              |
| JBH2245 | LiNX DLX-IDC (Intuitive Dual Control)                                                                               |
| JBH2259 | Dual Control Standard Bracket RH                                                                                    |
| JBH2260 | Dual Control Standard Bracket LH                                                                                    |
| JBH1221 | IDC Dual Control Detachable Bracket RH                                                                              |
| JBH1222 | IDC Dual Control Detachable Bracket LH                                                                              |
| JBH9260 | Bracket for mo-vis Scoot                                                                                            |
| JBH2231 | LiNX DLX-CR400 (Compact Remote)                                                                                     |
| JBH2232 | LiNX DLX-CR400LF (Compact Remote Low Force)                                                                         |
| JBH2207 | RH Mount                                                                                                            |

|         |                                                                                                                                                                                                                          |
|---------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| JBH2208 | LH Mount                                                                                                                                                                                                                 |
| JBH2265 | Centre Mount Swing Away to RH (Nucleus)                                                                                                                                                                                  |
| JBH2266 | Centre Mount Swing Away to LH (Nucleus)                                                                                                                                                                                  |
| JBH3410 | mo-vis Manual Lateral Swing Away RH Mounted                                                                                                                                                                              |
| JBH3411 | mo-vis Manual Lateral Swing Away LH Mounted                                                                                                                                                                              |
| JBH3412 | mo-vis Manual Vertical Swing Away RH Mounted                                                                                                                                                                             |
| JBH3413 | mo-vis Manual Vertical Swing Away LH Mounted                                                                                                                                                                             |
| JBH3420 | mo-vis Powered Lateral Swing Away RH Mounted                                                                                                                                                                             |
| JBH3421 | mo-vis Powered Lateral Swing Away LH Mounted                                                                                                                                                                             |
| JBH3422 | mo-vis Powered Vertical Swing Away RH Mounted                                                                                                                                                                            |
| JBH3423 | mo-vis Powered Vertical Swing Away LH Mounted                                                                                                                                                                            |
| JBH3200 | mo-vis Proportional Head Control Omni                                                                                                                                                                                    |
| JBH3221 | mo-vis Switched Head Control Omni                                                                                                                                                                                        |
| JBH3222 | mo-vis Switched Head Control w/ Sip & Puff Omni                                                                                                                                                                          |
| JBH2217 | Bracket RH, mounted for user                                                                                                                                                                                             |
| JBH2218 | Bracket LH, mounted for user                                                                                                                                                                                             |
| JBH2291 | Bracket RH, mounted for attendant                                                                                                                                                                                        |
| JBH2292 | Bracket LH, mounted for attendant                                                                                                                                                                                        |
| JBH2273 | DLX-IN200 Module (Jack port for external power button)                                                                                                                                                                   |
| JBH2254 | Egg Switch, black                                                                                                                                                                                                        |
| JBH1711 | Red Piko Button                                                                                                                                                                                                          |
| JBH1712 | Yellow Piko Button                                                                                                                                                                                                       |
| JBH1713 | Green Piko Button                                                                                                                                                                                                        |
| JBH1714 | Orange Piko Button                                                                                                                                                                                                       |
| JBH1715 | Blue Piko Button                                                                                                                                                                                                         |
| JBH1716 | Black Piko Button                                                                                                                                                                                                        |
| JBH4300 | Hand Warmer                                                                                                                                                                                                              |
| JBZ4301 | Hand Warmer 2 Dome Hand Hood                                                                                                                                                                                             |
| JBZ4302 | Hand Warmer 1 Dome Hand Hood                                                                                                                                                                                             |
| JBH4400 | Micro Joystick Ball                                                                                                                                                                                                      |
| JBH4401 | Micro Joystick Cup 10.1mm                                                                                                                                                                                                |
| JBH4402 | Micro Joystick Cup 15.5mm                                                                                                                                                                                                |
| JBH1294 | Indoor setting, reduced speed approx. 3 km/h                                                                                                                                                                             |
| JBH1908 | 10 km/h Package<br>Always combined with LED Light system with Indicators, 30 Lux headlights                                                                                                                              |
| JBH1909 | 12 km/h Ultra High Speed Package<br>Always combined with LED Light system with Indicators, 30 Lux headlights                                                                                                             |
| JBH1437 | LED Light system with indicators 30 Lux headlights, included with 10km/h motors                                                                                                                                          |
| JBH1254 | Without Light system                                                                                                                                                                                                     |
| JBH1244 | Kerb Rider with park position feature                                                                                                                                                                                    |
| JBH1450 | Preparation for DAHL-Docking station (Adequate headrest is compulsory).<br>Includes a mounting bracket in the battery box.<br>Docking Station for car implementation has to be ordered separately with DAHL Engineering. |
| JBH1401 | 6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19).<br>Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory.                           |
| JBH2894 | Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down.<br>For transport worthy versions (1494/1401) only                                                                             |

|         |                                                            |
|---------|------------------------------------------------------------|
| JBH1415 | Rear view mirror RH                                        |
| JBH1339 | Rear view mirror LH                                        |
| JBH2128 | Bag Hooks (Max load up to 6kg = impact on max user weight) |
| JBH1312 | Cane holder RH                                             |
| JBH1363 | Cane holder LH                                             |
| JBH1416 | Fender for Castors                                         |
| JBH1404 | Tray, Swing Away RH                                        |
| JBH1405 | Tray, Swing Away LH                                        |



For a more complex customer requirement,  
please email: [salesUK@dhcare.com](mailto:salesUK@dhcare.com)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):