



Yes, you can:

Tel: 0800 043 0881
Email: salesUK@dhcare.com
Online Spares available at: www.dhcare.co.uk

Power Retail Prescription Form

September 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Erivo

Standard Features

- Max user weight 150kg (23.5 Stone)
- 6kmh
- High quality foam seat cushion
- Solid seat plate
- Padded backrest with large built in pocket
- Fold-up armrests
- Fold up and height adjustable footplate (3 positions)
- LiNX DLX-REM060 remote with simple swivel arm
- 8" PU front castor
- Rear wheel suspension with 12.5" drive wheel, solid tyre
- 2 Year warranty



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL		
ERIR10/11	Erivo Stock	☐
ERIVO		
1595187	R10, 11Ah Battery, 17KM range	£2,900
1595188	R11, 30Ah battery, with lighting package, 35KM range	£3,500
ACCESSORIES (To be sent as spare parts for retrofitting. Additional lead times possible)		
ERI1234	1594974 - Bracket Kit for Attendant Mounting (includes 1000mm bus cable to fit REM 060 or REM 050 in the attendant position)	£196
ERI5678	1602792 - REM 050 Dual Control (to work in conjunction with REM 060 - Bracket Kit for Attendant 1594974 is required)	£192

ACCESSORIES SHIPMENT*

DNPS987		DO NOT PART SHIP - Accessories to be shipped with Erivo	◆
---------	--	---	---

	TOTAL PRICE	
--	-------------	--

Non Conforming Options:



For a more complex customer requirement,
please email: salesUK@dhcare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):