



Yes, you can:

Tel: 0800 043 0881
Email: salesUK@dhcare.com
Online Spares available at: www.dhcare.co.uk

Power Retail Prescription Form

September 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX10 Stock

Standard Features

- Max user weight 136kg (21.4 Stone)
- Crash Tested ISO-7176
- 6kph
- Ultra compact base width
- Highly cross-compatible
- Action style seating
- Solid seat plate
- 2 seat widths (430mm / 480mm) with 50mm width adjustment
- 410mm - 510mm seat depth adjustment (460mm standard)
- TAB Backrest with nylon overlay, 510mm high
- Height, depth, width adjustable Action style arm rests and pads
- Standard 80 Action style hangers, length adjustable with flip up footplates
- Drive Only & manual recline / Powered Tilt & manual recline / Powered Tilt & Powered Recline SKU's
- LINX DLX-REM060 or DLX-REM210 remote with R/H swing away bracket (bracket can be L/H by reconfiguring)
- 8" solid castors
- Rear wheel suspension with 14" solid drive wheels
- 51Ah GEL batteries
- 2 seat to ground heights (450mm standard / 470mm)
- Neckrest included with powered seating SKUs
- Charger & User Manual included



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code **Description**

BASE MODEL

AVIRX10	☐	AVIVA RX10 Stock	☐
---------	--------------------------	------------------	-----------------------

AVIVA RX10 Drive Only		
1680059		Drive Only / 430mm Seat Width / 460mm Seat Depth / 6kph / 450mm STGH / DLX-REM060 / RH Swing Away Bracket £3,650
1680060		Drive Only / 480mm Seat Width / 460mm Seat Depth / 6kph / 450mm STGH / DLX-REM060 / RH Swing Away Bracket £3,650
AVIVA RX10 Powered Tilt		
1680069		Powered Tilt / 430mm Seat Width / 460mm Seat Depth / 6kph / 450mm STGH / DLX-REM210 / RH Swing Away Bracket / Neckrest £3,800
1680070		Powered Tilt / 480mm Seat Width / 460mm Seat Depth / 6kph / 450mm STGH / DLX-REM210 / RH Swing Away Bracket / Neckrest £3,800
AVIVA RX10 Powered Tilt & Powered Recline		
1680210		Powered Tilt & Powered Recline / 430mm Seat Width / 460mm Seat Depth / 6kph / 450mm STGH / DLX-REM210 / RH Swing Away Bracket / Neckrest £3,950
1680211		Powered tilt & Powered Recline / 480mm Seat Width / 460mm Seat Depth / 6kph / 450mm STGH / DLX-REM210 / RH Swing Away Bracket / Neckrest £3,950

TOTAL PRICE	
-------------	--

Non Conforming Options:



For a more complex customer requirement,
please email: salesUK@dhcare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):